# lonizer

# **INSTRUCTION MANUAL**



### ©2021 A&D Company Ltd. All rights reserved.

- □ No part of this publication may be reproduced, transmitted, transcribed, or translated into any language in any form by any means without the written permission of A&D Company Ltd.
- □ The contents of this manual and the specifications of the instrument covered by this manual are subject to change for improvement without notice.



1WMPD4004479A

# Introduction

This manual describes outline of the AD-1683A ionizer and how to use it correctly. Please thoroughly read this manual before using the ionizer and keep it at hand for future reference.

#### Features 2.

This ionizer eliminate static electricity on the weighing sample by generating bipolar ions from 4 electrode needles.

- □ The ionizer can eliminate static electricity on a charged weighing sample without breeze. A weighing error can be reduced.
- □ This ionizer can be controlled by using built-in IR sensor, option IR switch (AX-IR-SWITCH).
- Electrode unit is removable. A unit can be cleaned and replaced.


### Static electricity

In general, nonconductors such as powder, filter, medicine wrapping paper, plastic etc. easily become electrostatically charged when the ambient humidity is less than 45 %RH. The static electricity may cause weighing error of approximately few mg at weighing. This ionizer can perform static elimination effectively.

### Cautions before Use 3.

- Read the following cautions for safety use of the ionizer.
- □ Operate the ionizer gently because of precision instrument.
- □ Do not install the following place.
- Place getting water, vibration, shock, direct sunshine, Dusty place, air including salt or corrosive gas, a place in flammable gas.
- Do not turn on power of the ionizer until the installation is finished. The switch to turn off is not equipped in the ionizer.

4. Specification	าร	
lon generation method	Direct corona discharge	
Effective range of static electricity elimination	Approximately 7 cm or shorter from the electrode needle	
Operation temperature and humidity	5 to 40 °C, 85% RH or less (no condensation)	
Discharge electrode needle	Tungsten (with a life of approx. 10,000 hours)	
Dimensions	68(W)×128(D)×163(H) mm	
Mass	Approximately 370 g	



### Caution

6.

6.2.

elimination starts.

Confirm that voltage, frequency and outlet type is correct for your local voltage.

6.1. Preparing the AC adapter

How to Use

1. Connect the accessory AC adapter to

2. When the AC adapter is connected to

the outlet, the power lamp lights up.

Static elimination

Example of the Outlet

AC adapter jack

1. Put the weighing sample into the effective range shown in the figure 1.

front panel, option IR switch, the ION lamp lights up and static

3. Perform static elimination referring to the figure 1. Static elimination


factory setting. Referring to "7. Optimizing the lonizer", optimize


electrode needles and weighing sample is 10 cm or father.

process stops automatically in 3 seconds and ION lamp turns off for

Example of the AC adapter

the AC adapter jack.






Caution Remove power supply cable when connecting peripherals to the ionizer.

Static elimination process can be controlled by using option IR switch instead of the built-in IR sensor. The ionizer cannot send PRINT and RE-ZERO command to the balance by using option IR switch.

Option IR switch

# 7. Optimizing the lonizer

Static elimination method can be optimized by using switches on the rear. Refer to the following table for functions of switches. Change the settings of the switch 1 and switch 2 eliminating static electricity outside the effective range shown as in the figure 1.



Factory settings of switches is all ON. Knob of switch

Factory settings of switches is all ON (Upper side). It means stat elimination method is "Timer mode", "Discharging time" is 3 second "Built-in IR sensor" and "Buzzer" is available

Table 1 Fu	unctions of the	optimizing	switches
------------	-----------------	------------	----------

Table 1 1 unctions of the optimizing switches						
Switch No. / Item		of switch	Descriptions			
Switch 1		Upper side	Timer mode #	1		
Static elimination method	OFF	Lower side	Manual mode #	2		
Switch 2	ON	Upper side	3 seconds			
Discharging time #3	OFF	Lower side	10 seconds			
Switch 3	ON	Upper side	Available			
Built-in IR sensor	OFF	Lower side	Not Available			
Switch 4	ON	Upper side	Available			
Buzzer #4	OFF	Lower side	Not Available			

When switch 1 is selected "Timer mode", static elimination is performed for "Discharging time" set by switch 2.

- Once static elimination has started, it continues until built-in IR 2. Responding the IR sensor (touchless infrared proximity sensor) on the #2 sensor or option IR switch responds again when switch 1 is selected "Manual mode". "Manual mode" does not stop static elimination automatically.
  - The setting of the switch 2 is effective when switch 1 is selected "Timer mode".
  - static elimination method and discharging time if the distance between #4 Buzzer sounds at turning the power on or responding the built-in IR sensor or option IR switch.




tic	
onds,	
	IF

#### Maintenance 8.

- **Caution** Do not touch the discharge electrode unit to avoid electric shock during neutralization.
- □ When the ionizer is used continuously, discharge electrode needle and around may get dirty and neutralization performance may become weak. Clean the electrode needle periodically using cotton swabs to maintain performance.
- □ When the electrode needle wears out and static elimination performance does not refresh by cleaning, replace all of discharge electrode units to new ones of option. Life time of discharge electrode unit is approximately 10000 hours.

## **Procedure of replacement**

- 1. Remove the AC adapter to turn off the ionizer.
- 2. Rotate a discharge electrode unit to 45 degree counterclockwise. Remove it.
- 3. Insert new units and rotate them to 45 degree clockwise.



Discharge electrode unit

# 9. Options

## Discharge electrode unit

- AX-BM-NEEDLESET (A set is 4 units.)
- □ Replacement electrode units.
- Replace all of 4 units with new ones.
- □ Refer to the "8. Maintenance" for replacement procedure.

## IR switch

AX-IR-SWITCH

- □ Touchless infrared proximity sensor.
- □ Refer to the "6.3. Connecting Options" how to connect it.
- □ Static elimination operation can be controlled by moving hand over the "SENSOR".



