


# Instruction Manual for AD-1675A Tabletop Breeze Break (S)

## 1. Features


- The tabletop breeze break protects the balance from breezes due to air conditioning and body motion, thereby reducing balance weighing error.
- The transparent panel assembly consists of antistatic plastic material, which protects the balance from static electricity.
- The hole on the back panel can be used to pass the RS-232C and AC adapter cables through.
- Sliding panels on top, right and left side enable samples to be inserted or removed from each side.

## 2. Specifications

- Weight : Approximately 2 kg
- Size : 257 (Width) x 371 (Depth) x 258 (Height) mm
- Material : Transparent panel : Antistatic plastic  
Frame : Aluminum
- Applicable models (balances made by A&D)
  - BA-T/BA series (1 $\mu$ g model)
  - AD-4212A series
  - AD-4212B series
  - AD-4212C series
  - AD-4212D series



## 3. External dimensions



### Caution

- Do not use organic solvents to clean the tabletop breeze break as that may damage the transparency and antistatic properties of the panels.
- Use a soft lint free cloth dampened with warm water and a mild detergent when cleaning the tabletop breeze break.