# FEEEE ESTRICT

# INSTRUCTION MANUAL

# **Carat Balance**

FZ-1200CT / FZ-700CT / FZ-500CT FX-1200CT / FX-700CT / FX-500CT



© 2012 A&D Company, Limited. All rights reserved.  No part of this publication may be reproduced, transmitted, transcribed, or translated into any language in any form by any means without the written permission of A&D Company, Limited.
The contents of this manual and the specifications of the instrument covered by this manual are subject to change for improvement without notice.
Windows, Word and Excel are the registered trademarks of the Microsoft Corporation.
The contents of this manual and the specifications of the instrument covered by this manual are subject to change for improvement without notice.

# Contents

	peration	
1.	Introduction	
1-1.	About This Manual	
1-2.	Features	
1-3.	Compliance	4
2.	Unpacking The Balance	<u>5</u>
 2-1.	Unpacking	
2-2.	Installing The Balance	
3.	Precautions	7
o. 3-1.	Before Use	
3-2.	During Use	
3-3.	After Use	
3-4.	Power Supply	
4.	Display Symbols And Key Operation	
5.	Weighing	
5-1.	Basic Operation (Gram Mode)	11
Adapti	ng To The Environment	
6.	Response Adjustment	12
7.	Calibration	13
7-1.	Calibration Mode	
7-2.	Calibration Using The Internal Mass (One-Touch Calibration, Only For The FZ-CT Series)	
7-3.	Calibration Test Using The Internal Mass (Only For The FZ-CT Series)	
7-4.	Calibration Using An External Weight	
7 <b>-</b> 5.	Calibration Test Using An External Weight	
7 <b>-</b> 6.	Correcting The Internal Mass Value 1 (Only For The FZ-CT Series)	
7-7.	Correcting The Internal Mass Value 2 (Only For The FZ-CT Series)	
F ati	· · · · · · · · · · · · · · · · · · ·	
Function		20
8.	Function Switch And Initialization	
8-1.	Permit Or Inhibit	
8-2.	Initializing The Balance	
9.	Function Table	
9-1.	Structure And Sequence Of The Function Table	
9-2.	Display And Keys	
9-3.	Details Of The Function Table	
9-4.	Description Of The Class "Environment, Display"	
9-5.	Description Of The Item "Data Output Mode"	
9-6.	Description Of The Item "Data Format"	
9-7.	Data Format Examples	
9-8.	Weighing Units	
9-9.	Storing Units	
9-10.	<b>5</b>	
9-11.		
9-12.	· · · · · · · · · · · · · · · · · · ·	
9-13.	Statistical Calculation Mode	40

9-14.	Comparator Function	44
9-15.	Clock And Calendar Function (Only For The <b>FZ-CT</b> Series)	45
10. 10-1.	ID Number And GLP Report	
10-2.		
11.	Underhook	51
12.	Programmable-Unit	52
13.	Password Function	53
13-1.	Using Password Function	53
13-2.		
13-3.	1 9	
13-4.	Missing Password	55
Interfa	ce And Communication	
14.	RS-232C Serial Interface	56
15.	Connection To Peripheral Equipment	57
15-1.	Connection To The AD-8121B Printer	57
15-2.	· ·	
15-3.	Using Windows Communication Tools (WinCT)	58
16.	Commands	59
16-1.		
16-2.	<b> </b>	
16-3.	Settings Related To RS-232C	61
Mainte	nance	
17.	Maintenance	62
18.	Troubleshooting	63
18-1.	Checking The Balance Performance And Environment	63
18-2.		
18-3.	1 7	
18-4.	Asking For Repair	65
19.	Options And Peripheral Equipments	66
20.	Specifications	69
21.	External Dimensions	71
22.	Terms/Index	72
22-1.	Terms	72
22-2.	Index	73

# 1. Introduction

This manual describes how the FZ-CT/FX-CT series balance works and how to get the most out of it in terms of performance. Read this manual thoroughly before using the balance and keep it at hand for future reference.

## 1-1. About This Manual

Basic Operation ...... Describes precautions on handling the balance, balance

construction and basic balance operation.

Adapting To The Environment..... Describes response adjustment and calibration.

Functions ...... Describes various functions of the balance.

Interface And Communication ..... Describes the interface which transmits data and controls the

balance.

Maintenance ...... Describes maintenance, error codes, troubleshooting,

specifications and options.

#### 1-2. Features

- Compact general-purpose balance, can be installed almost anywhere.
- Large breeze break, provided as a standard accessory for more accurate weighing.
- The breeze break can be removed with one touch, so cleaning is simple.
- Casing construction, strong protection against dust and moisture.
- Carat pans for weighing jewelry, provided as standard accessories.
- Multiple weighing units with most of the common units used around the world.
   The default settings of the unit of weight are the metric carat (ct) and grams (g). The unit can be switched by pressing the MODE key.
- Standard RS-232C serial interface to communicate with a computer and to output the Good Laboratory Practice (GLP) / Good Manufacturing Practice (GMP) compliant report.
- Statistical calculation mode to statistically calculate the weight data, and display or output the sum, maximum, minimum, range (maximum-minimum), average, standard deviation and coefficient of variation.
- Comparator Indicators, displaying the comparison results.
- Hold Function, provided for weighing a moving object such as an animal.
- As options, the USB interface ( $FX_{i}-02$ ), the Ethernet interface ( $FX_{i}-08$ ) and the built-in battery unit ( $FX_{i}-09$ ) are available.
  - Connecting the USB interface, by cable to a Windows computer, allows transmission of the weight data to Excel or Word. The Ethernet interface allows connection of the balance to a personal computer through a LAN. The battery unit allows use of the balance where a power source is not available.
- □ The AD-8920A remote display, that is sold separately, allows reading the weight value in a place away from the balance.

#### 1-3. Compliance

#### **Compliance with FCC Rules**

Please note that this device generates, uses and can radiate radio frequency energy. This device has been tested and has been found to comply with the limits of a Class A computing device pursuant to Subpart J of Part 15 of FCC rules. These rules are designed to provide reasonable protection against interference when this device is operated in a commercial environment. If this unit is operated in a residential area, it may cause some interference and under these circumstances the user would be required to take, at his own expense, whatever measures are necessary to eliminate the interference. (FCC = Federal Communications Commission in the U.S.A.)

#### Compliance With EMC Directives

This device features radio interference suppression and safety regulation in compliance with the following Council Directives

Council directive 2004/108/EC EN61326 EMC directive

Council directive 2006/95/EC EN60950 Safety of Information Technology Equipment

The CE mark is an official mandatory European marking.

Please note that any electronic product must comply with local laws and regulations when sold or used anywhere outside Europe.




# 2. Unpacking The Balance

# 2-1. Unpacking

- □ The balance is a precision instrument. Unpack the balance carefully. Keep the packing material to be used for transporting the balance in the future.
- See the illustrations to confirm that everything is included.

#### Install the large breeze break when drafts are present



# 2-2. Installing The Balance

Install the balance as follows:

- 1. Place the balance on a solid weighing table. Refer to "3. Precautions" for installing the balance.
- 2. Assemble the pan support, weighing pan, breeze break ring and breeze break, on the balance as shown in the illustration on page 5.
- 3. Adjust the leveling feet to level the balance. Confirm it using the bubble spirit level.
- 4. Confirm that the AC adapter type is correct for the local voltage and power receptacle type.
- 5. Connect the AC adapter to the balance. Warm up the balance for at least one hour with nothing on the weighing pan.

#### **Note**


The AC adapter plug is protected against dust and may be difficult to insert. When inserting the plug, turn the plug while pushing on it.


# 3. Precautions

To get the optimum performance from the balance and acquire accurate weight data, note the following:

#### 3–1. Before Use

- □ Install the balance in an environment where the temperature and humidity are not excessive. The best operating temperature is about 20°C / 68°F at about 50% relative humidity.
- Install the balance where it is not exposed to direct sunlight and it is not affected by heaters or air conditioners.
- Install the balance where it is free of dust.
- Install the balance away from equipment which produces magnetic fields.
- Install the balance in a stable place avoiding vibration and shock. Corners of rooms on the first floor are best, as they are less prone to vibration.
- □ The weighing table should be solid and free from vibration, drafts and as level as possible.
- Level the balance by adjusting the leveling feet and confirm it using the bubble spirit level.






- Ensure a stable power source when using the AC adapter.
- Plug in the AC adapter and warm up the balance for at least one hour.
- Calibrate the balance before use or after having moved it to another location.

Caution: Do not install the balance where flammable or corrosive gas is present.

# 3-2. During Use

- Discharge static electricity from the material to be weighed (hereinafter referred to as the sample). When a sample could have a static charge, the weight data is influenced. Ground the balance and try the following:
  - Eliminate static electricity using the optional AD-1683 DC static eliminator.
  - Try to keep the ambient humidity above 45%RH.
  - Use a metal shield case for a charged sample.
  - Wipe charged plastic samples with a moistened cloth.
- This balance uses a strong magnet as part of the balance assembly, so please use caution when weighing magnetic materials such as iron. If there is a problem, use the underhook on the bottom of the balance to suspend the material away from the influence of the magnet.
- Eliminate the temperature difference between a sample and the environment. When a sample is warmer (cooler) than the ambient temperature, the sample will be lighter (heavier) than the true weight. This error is due to a rising (falling) draft around the sample.
- Make each weighing gently and quickly to avoid errors caused by the evaporation of moisture from the sample or the absorption of moisture by the sample. Use a pair of long tweezers to place a sample on the pan to avoid placing a hand inside the breeze break.
- Do not drop things upon the weighing pan, or place a sample on the pan that is beyond the balance weighing capacity.
   Place a sample in the center of the weighing pan.
- Do not use a sharp instrument such as a pencil to press the keys. Use your finger only.
- Press the RE-ZERO key before each weighing to eliminate possible errors.
- Calibrate the balance periodically so as to eliminate possible errors.
- □ Take into consideration the affect of air buoyancy on a sample when more accuracy is required.
- Keep the balance interior free of dust and foreign materials.
- The breeze break (The FZ-CT /FX-CT series) is provided as a standard accessory. An anti-static treatment has been applied to the breeze break components, but they may be charged with static electricity when they are unpacked or when the humidity is low. If the weight value is unstable even when there are no drafts or the balance has a problem with repeatability, remove the breeze break. Or wipe the clear plates with a moistened cloth, use the optional AD-1683 DC static eliminator or apply an anti-static spray.



#### 3-3. After Use

- Avoid mechanical shock to the balance.
- □ Do not disassemble the balance. Contact the local A&D dealer if the balance needs service or repair.
- Avoid dust and water so that the balance weighs correctly. Protect the internal parts from liquid spills and excessive dust.
- Do not use organic solvents or chemical cloth to clean the balance.
   Clean the balance as described below:

· ·	Use a soft, lint free cloth that is moistened with a mild detergent to clean.  The edge of the weighing pan is sharp. Use caution when cleaning the pan.
Breeze break (standard accessory)	An antistatic treatment has been applied to the breeze break components. Use a soft, dry, lint free cloth to clean. Cleaning with a cloth that is moistened with water or a mild detergent repetitively, or washing with water, may degrade the antistatic effect.

# 3-4. Power Supply

When the AC adapter is connected, the balance is in the standby mode if the standby indicator is on (refer to "4. Display Symbols And Key Operation"). This is a normal state and does not harm the balance. For accurate weighing, allow the balance to reach equilibrium with room temperature, and then warm up the balance for at least one hour before use.

# 4. Display Symbols And Key Operation

#### **Key Operation**


Key operation affects how the balance functions. The basic key operations are:

- "Press and release the key immediately" or "Press the key" ......
  - = normal key operation during measurement
- "Press and hold the key"



#### **Display Symbols**

- Displays the weight data relative to the weighing capacity, in percentage, in the weighing mode (Capacity indicator)
- Number of statistical data (Statistical calculation mode)
- Displays the setting value in the function mode



Each key, when pressed or when pressed and held, functions as follows:

Key	When pressed	When pressed and held	
I/O ON:OFF	□ Turns the display on or off. The standby indicator is displayed when the display is turned off. The weighing mode is enabled when the display is turned on. When the display is turned on, the password input ( PR55 □ □ ) is displayed if the password function is used. The weighing mode is enabled if the password function is not used. Refer to "13–3. Inputting Password When Turning On The Balance".  □ This key is available anytime. Pressing the key during operation will interrupt the operation and turn the display off.		
1/10d SAMPLE	<ul> <li>In the weighing mode, turns the minimum weighing value on or off.</li> <li>In the counting or percent mode, enters the sample storing mode.</li> </ul>	Enters the function table mode. Refer to "9. Function Table".	
MODE	Switches the weighing units stored in the function table. Refer to "9-8. Weighing Units".	Enters the response adjustment mode.	
CAL	<ul> <li>Cancels the operation when performing settings. With the FZ-CT series, enters the calibration mode using the internal mass. (One-touch calibration)</li> </ul>	Enters the calibration mode.	
Q PRINT	<ul> <li>Outputs the weight data to a printer or personal computer using the RS-232C serial interface, depending on the function table settings.</li> <li>Confirms the operation when performing settings.</li> </ul>	<ul> <li>No function at the factory setting</li> <li>By changing the function table:</li> <li>Outputs "Title block" and "End block" for GLP/GMP compliant report. (Refer to "10-2. GLP Report".)</li> </ul>	
→0← RE-ZERO	Sets the display to zero.		

# 5. Weighing

The default settings of the unit of weight are the metric carat (z t) and grams (g).

Pressing the MODE key switches the displayed unit, between £t and 9.

To use other units, refer to "9–9. Storing Units" on page 37 to store the units.

## 5–1. **Basic Operation** (Gram Mode)


- 1 Place a container (Carat pan) on the weighing pan, if necessary.


  Press the RE-ZERO key to cancel the weight (tare). The balance displays QQQQ =t.
- 2 Place a sample on the pan or in the container.
- Wait for the stabilization indicator to turn on. Read the value. While the stabilization indicator is on, pressing the PRINT key will output the weight value, using the RS-232C serial interface.
- 4 Remove the sample and container from the pan.

#### Note

- Peripheral equipment, that is sold separately, such as a printer or a personal computer is required.
- When other units of weight are stored as described in "9–9. Storing Units", pressing the MODE key switches the displayed unit.
- Press the SAMPLE key to turn on or off the minimum weighing value.

When the ON:OFF key is pressed with a container placed on the weighing pan, the balance displays Outline and weighing is started,





# 6. Response Adjustment

This function stabilizes the weight value, reducing the influence on weighing that is caused by drafts and vibration at the place where the balance is installed.

The function has three stages as follows and can be changed by simple key operation.

Indicator	Description	
FAST	Fast response, but prone to drafts and vibration.	
FASI	Good for target weighing.	
MID.	•	
SLOW	Slow response, but strong against drafts and vibration	
SLOW	Good for weighing which requires a stable display.	



#### Operation


- 1 Press and hold the MODE key until RESPONSE is displayed. And then, release the key.
- 2 Press the MODE key to select a rate of the response adjustment. Either FAST, MID, or SLOW can be selected.
- 3 After a few seconds of inactivity, or when the PRINT key is pressed, the balance displays End. Then, it returns to the weighing mode and displays the updated response indicator. The response indicator remains displayed for a while.

#### Note

□ Setting the response rate automatically changes the "Condition ([and)" and "Display refresh rate (5Pd)" parameters of "Environment, Display (bR5Fnc)" in the function table, as shown below:

Indicator	Cond	(Condition)	5Pd (Display refresh rate)	
FAST		0	(Approx. 10 times/second)	
MID.		1	[] (Approx. 5 times/second)	
SLOW		2	☐ (Approx. 5 times/second)	

When the balance is to be used with other setting combinations, set each parameter in the function table.



# 7. Calibration

## 7-1. Calibration Mode

The balance has the following modes. (Calibration test does not perform calibration.)

- Calibration using the internal mass (One-touch calibration, only for the FZ-CT series)
- Calibration using an external weight
- Calibration test using an internal mass (Only for the FZ-CT series)
- Calibration test using an external weight

#### Terms

The following terms are defined as follows:

External weight = A weight that you have. Referred to as a calibration weight when used for

calibration.

Calibration weight = A weight used for calibration

Target weight = An external weight used for calibration test

#### Caution

When calibrating, be sure to install the provided large breeze break.

Calibration adjusts the balance for accurate weighing.

Besides periodic calibration and before each use, perform calibration when:

- The balance is installed for the first time.
- The balance has been moved.
- The ambient environment has changed.
- Do not allow vibration or drafts to affect the balance during calibration.
- To output the GLP/GMP compliant report using the RS-232C serial interface, set "GLP output ( μη Γα )" of "Data output ( dαμ \( ) \)". For details, refer to "9. Function Table". To display the GLP output, a printer or personal computer is required.
- □ Calibration test is available only when "GLP output ( ייה + □ )" of "Data output ( dout )" is set to " /" or "?".

## Caution On Using An External Weight

The accuracy of an external weight will influence the accuracy of weighing. Select an appropriate weight as listed below:

Model	Usable Calibration Weight	Adjustable Range
FZ-1200CT / FX-1200CT	250 g, <b>200 g</b> , 100 g, 50 g, 20 g	
FZ-700CT/ FX-700CT	150 g, <b>100 g</b> , 50 g, 20 g	-0.0150 g to +0.0150 g
FZ-500CT/ FX-500CT	<b>100 g</b> , 50 g, 20 g	

The calibration weight in bold type: factory setting.

The calibration weight value can be adjusted within the range above.

## Display

This indicator means "the balance is measuring calibration data".

Do not allow vibration or drafts to affect the balance while this indicator is displayed.

## 7–2. Calibration Using The Internal Mass (One-Touch Calibration, Only For The FZ-CT Series)

This function calibrates the balance using the internal mass. The only operation required is to press the CAL key. (The FX-CT series does not have the one-touch calibration function.)

#### Note

- Level the balance by adjusting the leveling feet and confirm that the bubble is in the center of the spirit level.
- Poor level adjustment may cause a calibration error while using the internal mass.

#### Operation

- 1 Connect the AC adapter and warm up the balance for at least one hour with nothing on the weighing pan.
- 2 Press the CAL key.
- 3 The balance displays EBL in and performs calibration using the internal mass. Do not allow vibration or drafts to affect the balance.
- 4 The balance displays  $\[ \underline{\mathcal{E}}_{nd} \]$  after calibration. If the "GLP output  $(\[ \underline{\mathit{InF}}_{0} \])$ " parameter of the function table is set to " $\[ \underline{\mathit{I}} \]$ ", the balance displays  $\[ \underline{\mathit{ILP}} \]$  and outputs "Calibration Report" using the RS-232C interface. For details on the calibration report format, refer to "10–2. GLP Report".
- 5 The balance will automatically return to the weighing mode after calibration.

#### **About The Internal Mass**

The mass of the internal mass may change due to corrosion or other damage caused by the operating environment, or due to aging. The internal mass value can be corrected using an external mass. To maintain the weighing accuracy, perform the internal mass value correction periodically. For details, refer to "7–4. Correcting The Internal Mass Value 1" and "7–7. Correcting The Internal Mass Value 2".

# 7–3. Calibration Test Using The Internal Mass (Only For The FZ-CT Series)

- This function tests the balance accuracy using the internal mass.
- Calibration test does not perform calibration.
- □ This is available only when the "GLP output ( ¬¬F□)" parameter is set to " /" or "²".


#### Operation

- 1 Connect the AC adapter and warm up the balance at least one hour.
- 2 Refer to "9. Function Tablel" to set "GLP output ( nFa)" to " l" or "l"
- 3 Press and hold the CAL key until [[[ in ] is displayed, then release the key.
- 4 The balance measures the zero point.

  Prevent vibration and drafts to affect the balance.
- 5 The measured zero point data is displayed.
- 6 The balance measures the internal mass.Prevent vibration and drafts to affect the balance.
- 7 The value of the internal mass is displayed. The normal range of the value is as follows:

Model	The Internal Mass	The Normal Range
FZ-1200CT	200.000 g	
FZ-700CT	100 000 a	±0.002 g
FZ-500CT	100.000 g	

- 8 If GLP output is set, ΔLP is displayed, the calibration test report is output to the RS-232C interfac. Refer to "GLP output ( ισΕα)" of the function table, "10-2. GLP Report". Εσd is displayed after the calibration.
- 9 The balance automatically returns to the weighing mode.



0000 st

## 7-4. Calibration Using An External Weight

This function calibrates the balance using an external weight.

#### Operation

- 1 Connect the AC adapter and warm up the balance for at least one hour with nothing on the pan.
- 2 Press and hold the CAL key until [FIL out] is displayed, and then release the key.
- 3 The balance displays [FIL []].
  - If you want to change the calibration weight (a list of usable weights is shown on page 13), press the SAMPLE key and proceed to step 4.
  - If you use the calibration weight value stored in the balance, proceed to step 5.
- 4 Specify the calibration weight value as follows:

SAMPLE key To switch the display condition to: "All of the segments blinking" (calibration weight selection mode) or "The last three digits blinking" (value adjustment mode).

RE-ZERO key To increase the value.

MODE key To decrease the value.

To select the calibration weight or adjust

the value.


PRINT key To store the new weight value. Even if the

AC adapter is removed, the data is maintained in non-volatile memory.

CAL key To cancel the operation and return to [FRL []].

5 Confirm that there is nothing on the pan and press the PRINT key. The balance measures the zero point. Do not allow vibration or drafts to affect the balance. The balance displays the calibration weight value.

- 6 Place the displayed calibration weight on the pan and press the PRINT key. The balance measures the calibration weight. Do not allow vibration or drafts to affect the balance.
- 7 The balance displays  $\boxed{\textit{End}}$ . Remove the weight from the pan.
- 8 When the GLP output is set, the balance displays LLP and outputs "Calibration Report".
- 9 The balance automatically returns to the weighing mode.
- 10 Press the MODE key to select the gram mode. Place the calibration weight on the pan and confirm that the value displayed is within ±0.002 g of the calibration weight value. If it is not within the range, check the ambient conditions such as drafts and vibration, also check the weighing pan. Then, repeat steps 1 to 9.



## 7-5. Calibration Test Using An External Weight

This function tests the balance weighing accuracy using an external weight and outputs the result.

This is available only when the "GLP output (info)" parameter is set to "1" or "2". (Calibration test does not perform calibration.).

#### Operation

- 1 Connect the AC adapter and warm up the balance for at least one hour with nothing on the pan.
- 2 Press and hold the CAL key until [[ aut is displayed, and then release the key.
- 3 The balance displays [[[ ]] .
  - If you want to change the calibration weight (a list of usable weights is shown on page 13), press the SAMPLE key and proceed to step 4.
  - If you use the calibration weight value stored in the balance, proceed to step 5.
- 4 Specify the calibration weight value as follows:

SAMPLE key To switch the display condition to: "All of the segments blinking" (target weight selection mode) or "The last three digits blinking" (value adjustment mode).

RE-ZERO key To increase the value.

MODE key To decrease the value.

To select the target weight or adjust the


value.

PRINT key To store the new weight value. Even if the

AC adapter is removed, the data is maintained in non-volatile memory.

CAL key To cancel the operation and return to [[ [ ] ].

- 5 Confirm that there is nothing on the pan and press the PRINT key. The balance measures the zero point and displays the measured value. Do not allow vibration or drafts to affect the balance. The balance displays the target weight value.
- 6 Place the displayed target weight on the pan and press the PRINT key. The balance measures the target weight and displays the measured value. Do not allow vibration or drafts to affect the balance.
- 7 The balance displays *End*. Remove the weight from the pan.
- 8 When the GLP output is set, the balance displays LLP and outputs "Calibration Test Report".
- 9 The balance automatically returns to the weighing mode.




# 7-6. Correcting The Internal Mass Value 1 (Only For The FZ-CT Series)

- The FZ-CT series balance can correct the internal mass value using an external weight.
   Method 1: Calibrate the balance using the internal mass. Measure the external weight.
   Calculate the correction value. And store it in the balance.
- The adjustable range of the internal mass value is as follows:

Model	Correction Reference Value	Adjustable Range
FZ-1200CT	200 g	
FZ-700CT	100 a	-0.0150 g to +0.0150 g
FZ-500CT	100 g	

#### Operation

Example : Correcting the weight value by +0.0030 g in 200 g using the FZ-1200CT. If correcting the weight value by +0.0030 g in 100 g, use the correction value of +0.0060 g as



- 1 Perform calibration using the internal mass (one-touch calibration).Place the external weight on the weighing pan to obtain the correction value.
- 2 With the factory setting, the balance can not correct the internal mass value. Refer to "8. Function Switch And Initialization" and set the internal mass value correction switch to 1.
- In the weighing mode, press and hold the SAMPLE key to display  $6855 \, \text{nc}$ .
- Press the SAMPLE key several times until [5 in ] is displayed. If [5 in ] is not displayed, perform step 2.
- 5 Press the PRINT key.

Correct the internal mass value using the following keys:

RE-ZERO key... To increase the value by one.


MODE key..... To decrease the value by one.

PRINT key...... To store the correction value and display the next menu item of the

function table.

CAL key ..... To cancel the operation and return to the weighing mode.

- 6 Press the CAL key to return to the weighing mode.
- 7 Press the CAL key to perform calibration using the internal mass.
- 8 Confirm that the internal mass value has been corrected properly by placing the external weight on the weighing pan.
  If the internal mass value has not been corrected properly, repeat the procedure to adjust the correction value.



0,000 g

## 7–7. **Correcting The Internal Mass Value 2** (Only For The FZ-CT Series)

- The FZ-CT series balance can correct the internal mass value using an external weight. Method 2: Calibrate the balance using the external weight. The balance performs automatic calibration using the internal mass, corrects the internal mass value and stores it in the balance. The corrected mass value is maintained in non-volatile memory even if the AC adapter is removed.
- □ The usable calibration weights and the adjustable range are as follows:

Model	Correction Reference Value	Adjustable Range
FZ-1200CT	250 g, <b>200 g</b> , 100 g, 50 g, 20 g	
FZ-700CT	150 g, <b>100 g</b> , 50 g, 20 g	-0.0150 g to +0.0150 g
FZ-500CT	<b>100 g</b> , 50 g, 20 g	


The calibration weight in bold type: factory setting.

The calibration weight value can be adjusted within the range above.

#### **Operation**

- 1 With the factory setting, the balance can not correct the internal mass value. Refer to "8. Function Switch And Initialization" and set the internal mass value correction switch to 1.
- 2 In the weighing mode, press and hold the SAMPLE key to display  $685 \, \text{Fnc}$ .
- 3 Press the SAMPLE key several times until 5 in 2 is displayed. If 5 in 2 is not displayed, perform step 1.
- 4 Press the PRINT key to display [FRL []].
  Refer to "7-4. Calibration Using an External Weight", to perform the calibration.
- 5 After calibration, remove the weight.
  The balance displays [FERBY]. Press the PRINT key.
- 6 [FIL 557] is displayed and the balance corrects the internal mass value automatically
- 7 After correcting the internal mass value, the balance displays [FRL 17] and performs the calibration automatically using the corrected internal mass value.
- 8 The balance displays *ξnd* and returns to the weighing mode.
- 9 Confirm that the internal mass value has been corrected properly by placing the external weight used for the correction procedure on the weighing pan. If the internal mass value has not been corrected

If the internal mass value has not been corrected properly, repeat the procedure to adjust the correction value.




# 8. Function Switch And Initialization

## 8-1. Permit Or Inhibit


The balance stores parameters that must not be changed unintentionally. There are two switches for the purpose of protecting these parameters. Each switch can select either "permit" or "inhibit". The "inhibit" protects parameters against unintentional operations.

	FZ-CT series	FX-CT series
Function switch	- Calibration using the internal mass	<ul><li>Function table</li><li>Calibration using an external weight</li></ul>

#### Switches (FZ-CT Series)



#### Switches (FX-CT Series)



#### Operation

- 1 Press the ON:OFF key to turn off the display.
- 2 While pressing and holding the PRINT key and the SAMPLE key, press the ON:OFF key. The balance displays 75.
- 3 Press the PRINT key. Then the balance displays the function switches.
- 4 Set the switches using the following keys.

SAMPLE key .... To select the switch to change the parameter. The selected switch blinks.

RE-ZERO key ··· To change the parameter of the switch selected.

☐ : To inhibit changes. /: To permit changes

PRINT key ...... To store the new parameter and return to the weighing mode.

CAL key .......... To cancel the operation (  $\lfloor \mathcal{L}_{\Gamma} \rfloor$  is displayed.)

To return to the weighing mode, press the CAL key once again.

## 8-2. Initializing The Balance

This function returns the following parameters to factory settings.


- Calibration data
- Function table
- □ The sample unit mass value (counting mode), 100% reference mass value (percent mode)
- External calibration weight
- Function switch settings
- Statistical data
- Internal mass value (Only for the FZ-CT series)

#### Note

Be sure to calibrate the balance after initialization.

## Operation

- 1 Press the ON:OFF key to turn off the display.
  While pressing and holding the PRINT key and the SAMPLE key, press the ON:OFF key. The balance displays P5.
- 2 Press the SAMPLE key to display [[Lr]].
- 3 Press the PRINT key.
  To cancel this operation, press the CAL key.
- 4 Press the RE-ZERO key.
- 5 Press the PRINT key to initialize the balance.
  The balance will automatically return to the weighing mode.



# 9. Function Table


The function table reads or rewrites the parameters that are stored in the balance. These parameters stored, even if the AC adapter is removed, are maintained in non-volatile memory.

# 9-1. Structure And Sequence Of The Function Table

The function table menu consists of two layers. The first layer is the "Class" and the second layer is the "Item". Each item stores a parameter.

#### **Example**

This example sets "Auto print mode A" for "Data output mode" and "10 digits" for "Auto print difference".



## 9-2. Display And Keys

Display/Key	Description
•	The symbol o indicates that the parameter displayed is in effect.
1/10d SAMPLE	When pressed and held in the weighing mode, enters the function table mode. Selects the class or item in the function table mode.
+0+ RE-ZERO	Changes the parameter.
Q PRINT	When a class is displayed, moves to an item in the class. When an item is displayed, stores the new parameter and displays the next class.
CAL	When an item is displayed, cancels the new parameter and displays the next class.  When a class is displayed, exits the function table mode and returns to the weighing mode.

# 9-3. Details Of The Function Table

Class	Item and Parame	ter		Des	cription					
	г ,	0	Fast response, sensitive	value	FAST	\\(\lambda \)\(\lambda \)\(\la				
	Cand Condition	- /	1		MID.	With "Hold I", sets the averaging time.				
	Contaition	5	Slow response, stable va	alue	SLOW	averaging time.				
	St-6		Stable range is ±1 digit	The sta	abilization	indicator illuminates when				
	Stability band width	- /	1	the display fluctuation is within the range.						
	Cability Baria Widti	2	Stable range is ±3 digits	With "	daLd i",S€	ets the stabilization range.				
	HoL d	<b>-</b> []	OFF	/ whe <u>n stable</u> in animal						
	Hold function	1	ON	mode.	With "HoL	d /", HOLD turns on.				
		0	OFF							
	trc	- /	Normal	Keens	zero displ	av by tracking zero drift				
וחכר	Zero tracking	2	Strong	пссро	Keeps zero display by tracking zero					
<i>bR5Fnc</i> Environment		3	Very strong							
Display	5Pd	<b>-</b> []	Approx. 5 times/second		Period to r	efresh the display				
. ,	Display refresh rate	1	Approx. 10 times/second	<u> </u>	i choa to i	circon tric display				
	Pnt	<b>-</b> []	Point (.)	Decim:	rmat					
	Decimal point	1	Comma (,)	Dooili	imal point format					
	P-on	<b>-</b> []	OFF	Turns	urns on the weighing mode display					
	Auto display-ON	1	ON	the AC	C adapter is connected.					
	Poff	<b>-</b> []	OFF			play after 10 minutes of				
	Auto display-OFF	1	ON (10 minutes)	inactivi	ity.					
		<b>-</b> []	Displays	r not to display the						
	Display at start	1	Does not display	minimu	nimum weighing value at weighing s					
	PEE6	0	Does not sound		or not to sound the beep					
	Веер	- /	Sounds	on keys.						
[L RdJ * Clock			r to "9–10. Clock and ndar function".	Confirms and sets the time and date. The time and date are added to output data.						
Clock		<b>-</b> []	No comparison	uirio ai	ia date are	c added to output data.				
		- 1	Comparison, excluding "n	ear zer	n" when s	table value or overloaded				
	CP .	2	Comparison, including "ne							
	Comp mode	3	Continuous comparison							
		4	Continuous comparison							
[P Fnc	bEP_	<b>-</b> []	OFF		9	20.0				
Comparator	LO buzzer	- 1	ON							
	ьер-	<u> </u>	OFF							
	OK buzzer	- 1	ON							
	bEP-	<b>-</b> []	OFF							
	HI buzzer	-	ON							
CP H I										
Upper limit		Dofo	rto "O 14 Commorator	Euss	tion"					
EP Lo		reie	r to " <b>9-14. Comparator</b>	runct	LIUII .					
Lower limit										

• : Factory setting Digit is a unit of minimum weighing value

\* : Only for the FZ-CT series

Class	Item and Parame	ter		De	escription
		<b>-</b> []	Key mode		Accepts the PRINT key only when the display is stable.
		1	Auto print mode A (Reference = zero)		Outputs data when the display is stable and conditions of
		2	Auto print mode B (Reference = last st	able value	##P-P, ##P-b and the reference value are met.
	Pr Ł Data output mode	3	Stream mode	data at the specified display refresh : Used when connected to the A remote display.	
		4	Key mode B (Immediately)	Accepts to	the PRINT key regardless of the ondition.
		5	Key mode C (When stable)	the displa	the PRINT key immediately when ay is stable, or waits for the display ole when not.
		5	Interval output mode	е	Uses interval output mode.
	RP-P	<b>-</b> []	Plus only		Displayed value>Reference
	Auto print polarity	1	Minus only		Displayed value <reference< td=""></reference<>
	rate print polarity	2	Both		Regardless of displayed value
	RP-6	0	10 digits		Difference between reference
	Auto print	<b>-</b> /	100 digits		value and displayed value
	difference	2	1000 digits		value alla alepiayea value
dout		0	Every display refres		
		<b>-</b>	2 seconds		
Data output		2	5 seconds		
	Int	3	10 seconds		Interval time for the interval output
	Interval time	4	30 seconds		mode
		5	1 minute		(With <sup>₽</sup> r Ł ᢃ)
		5	2 minute		
		7	5 minute		
		8	10 minute		
		<b>-</b> []	No output		
	5-Ed	- 1	Time only		
	Time/Date output	3	Date only		
		4	Time and date		
	5- 1d	<b>.</b> []	No output		Selects whether or not the ID
	ID number output	1	Output		number is output.
	PUSE	<b>-</b> []	No pause	`	Selects the data output interval.
	Data output pause	i	Pause (1.6 seconds	5)	·
	유는 - F Auto feed	<b>-</b> []	Not used		Selects whether or not automatic
	AUTO 1660	i ■ []	Used		feed is performed.
	ınFa	<b>-</b> []	No output		Selects the output format for the
	GLP output	i 2	AD-8121 format	4	GLP/GMP compliant report.
	0 1		General data forma	ı	
	유 d Zero after output	• []	Not used		Adjusts zero automatically after data is output
	Zeio aitei output	İ	Used		of minimum woighing value

Factory setting Digit is a unit of minimum weighing value
 Only for the FZ-CT series

Class	Item and Parame	eter		De	escription					
			600 bps							
		1	1200 bps							
	6PS	• 2	2400 bps							
	Baud rate	3	4800 bps							
		4	9600 bps							
		5	19200 bps							
	btPr	<b>-</b> []	7 bits, even							
	Data bit, parity bit	- 1	7 bits, odd							
			8 bits, none							
5 iF		<b>-</b> []	CR LF		CR: ASCII code 0Dh					
Serial	Terminator		CR		LF: ASCII code 0Ah					
interface		<b>-</b> []	A&D standard format							
	E Y P E	-	DP format							
			KF format		Refer to "9-6. Description of the Item "Data Format"".					
	Data format	3	MT format		item "Data Format".					
			NU format		_					
			CSV format							
	E-UP	0	No limit		Selects the wait time to receive a					
	Timeout	• /	1 second		command.					
	Er[d	<b>-</b> []	No output		AK:ASCII code 06h					
-, ,	AK, Error code	i	Output	Α.						
	ble-unit (Multi-unit)	Sets	an arbitrary coefficient.	ailable only when programmable- it mode is selected.						
ปก เะ Unit	Սո ւէ Unit		r to "9-8 <b>. Weighing Uni</b>	ts".						
ਾਰ ID number		Refer to "10. ID Number And GLP Report"								

• : Factory setting Digit is a unit of minimum weighing value

\* : Only for the FZ-CT series

## Caution

- The balance may not transmit the data completely at the specified refresh rate, depending on the baud rate.
- □ When a refresh rate of 20 times/second is used, set the baud rate to 4800 bps or higher.

Class	Item and Parame	ter		Description							
	RPF	<b>-</b> []	Normal weighing mode								
	Application function	1	Capacity indicator								
	7 ipplication ranductr	2	Statistical calculation mode	tatistical calculation mode							
		<b>-</b> []	Number of data, sum								
RP Fnc Application	SERF	1	Number of data, sum, maxi (maximum-minimum)	mum, minimum, average, range							
	Statistical function mode output items	2		Number of data, sum, maximum, minimum, average, range maximum-minimum), standard deviation, coefficient of variation							
	, , , , , , , , , , , , , , , , , , ,	3	Number of data, sum, maximum, minimum, average, range (maximum-minimum), standard deviation, coefficient of varia relative error								
LacFnc		<b>-</b> []		on the balance, the password input							
Password fur	nction	1	Used (PR55 #3) is displayed.								
Lac na Password			Refer to "13-2. Changing	Password".							
[5 in 2 *	value correction: Meth		Refer to " <b>7. Calibration</b> ".	Displayed only when the function switch (internal mass value correction switch) is set to !							
Internal mass	value correction: Meth	od 2		Series and the series of							

• : Factory setting Digit is a unit of minimum weighing value

\* : Only for the FZ-CT series

## 9-4. Description Of The Class "Environment, Display"

#### Condition ([and)

This parameter is for sensitive response to the fluctuation of a weight value. Used for powder target weighing, weighing a very light sample or when quick response weighing is required.

After setting, the balance displays FAST.

This parameter is for stable weighing with slow response. Used to prevent a weight value.

This parameter is for stable weighing with slow response. Used to prevent a weight value from drifting due to vibration or drafts. After setting, the balance displays SLOW.

#### Note

□ With "Hold function ( $H_0Ld$ )" set to "ON ( I)", this item is used to set the averaging time.

#### Stability Band Width (5t-b)

This item controls the width to regard a weight value as a stable value. When the fluctuation per second is less than the parameter, the balance displays the stabilization indicator and outputs the data. The parameter influences the "Auto print mode".

5t-b  $\Box$  This parameter is for sensitive response of the stabilization indicator. Used for exact weighing.

5 $\dot{\epsilon}$  This parameter ignores slight fluctuation of a weight value. Used to prevent a weight value from drifting due to vibration or drafts.

#### Note

□ With "Hold function ( $H_0Ld$ )" set to "ON ( I)", this item is used to set the stabilization range.

#### Hold Function (Hold) (Animal Weighing Mode)

This function is used to weigh a moving object such as an animal. When the weight data is 0.0200 g or greater and the display fluctuation is within the stabilization range for a fixed period of averaging time, the processing indicator illuminates and the balance displays the average weight of the animal. When the animal is removed from the weighing pan, the display returns to zero automatically. This function is available only when the hold function parameter is set to " /" (the animal mode indicator <code>HOLD</code> illuminates) and any weighing unit other than the counting mode is selected.

The stabilization range and averaging time are set in "Condition (Eand)" and "Stability band width (5E-b)".

Weighing Range
0.0200 g or greater

	Averaging	Time												
Cond O	Eand (1) 2 seconds Faster													
Cond I	4 seconds	<b>\$</b>												
Cond 2	8 seconds	More accurate												

Stabi	ization R	ange
St-6 0	Lesser	6.25%
5t-b l	1	12.5%
St-6 2	Greater	16.7%

#### **Zero Tracking** (t r c)

This function tracks zero point drift caused by changes in the environment and stabilizes the zero point. When the weight data is only a few digits, turn the function off for accurate weighing.

#### Note

- Digit indicates a unit of minimum weighing value.
- Erc Π The tracking function is not used. Used for weighing a very light sample.
- The tracking function is used. Normal zero tracking.
- *Erc ≥* The tracking function is used. Strong zero tracking.
- *Erc ∃* The tracking function is used. Very strong zero tracking.

#### Display Refresh Rate (5Pd)

Period to refresh the display. This parameter influences "Stream mode" along with "Baud rate" and "Data output pause".

#### Note

This item is selected automatically when the response rate is changed.

#### Decimal Point (Pnt)

The decimal point format can be selected.

## Auto Display-On (P-an)

When the AC adapter is connected, the display is automatically turned on in the weighing mode, without the ON:OFF key operation. Used when the balance is built into an automated system. One hour warm up is necessary for accurate weighing.

## Auto Display-Off (PoFF)

This function automatically turns off the balance display, if the weighing value is 0.5 g or less and the balance is in the no key operation or no command has been received for approx. 10 minutes. (The balance has the same power requirement in the weighing state as in the off state {standby mode}.)

## 9-5. Description Of The Item "Data Output Mode"

The parameter setting of the "Data output mode ( $P_{\Gamma} L$ )" applies to the performance when the data is transmitted using the RS-232C serial interface.

#### **Key Mode**

When the PRINT key is pressed with the stabilization indictor turned on, the balance outputs the weight data and the display blinks one time.

Required setting doub Prb 0 Key mode

#### Auto Print Modes A And B

When the displayed value is stable and the conditions of "Auto print polarity", "Auto print difference" and reference value are met, the balance outputs the weight data.

When the PRINT key is pressed with the stabilization indictor turned on, the balance outputs the data and the display blinks one time.

Mode A:

Example "For outputting the weight value each time a sample is added (or removed), with

"Ar -d" set to " /" (to adjust zero after the data is output)."

Required setting doub Prt Auto print mode A (reference = zero)

daut RP-P Auto print polarity

doub Re-d | Adjusts zero automatically.

Mode B:

Example "For outputting the weight value while a sample is added."

Required setting doub Prt 2 Auto print mode B (reference = last stable value)

daut RP-P Auto print polarity

#### Stream Mode

The balance outputs the weight data at the specified display refresh rate, regardless of the display condition. The display does not blink in this mode.

Example "For monitoring data on the AD-8920A remote display

Required setting doub Prb 3 Stream mode

basence 5Pd Display refresh rate

5 F 625 Baud rate

#### Caution

The balance may not transmit the data completely at the specified display refresh rate and baud rate.
 Set the baud rate higher.

#### **Key Mode B**

When the PRINT key is pressed, the balance outputs the weight data, regardless of the display condition. The display does not blink in this mode.

Required setting doub Prt 4 Key mode B

#### **Key Mode C**

When the | PRINT | key is pressed with the stabilization indictor turned on, the balance outputs the weight data. When the stabilization indicator is not on, the balance waits for the indicator to turn on and outputs the data. The balance blinks one time in this mode.

Key mode C

Prt 5 Required setting dout

#### Interval Output Mode

The weight data is periodically output. When the | PRINT | key is pressed, the balance starts to output the weight data at a preset interval time. When the | PRINT | key is pressed again, the balance stops outputting the weight data.

"For outputting the weight data periodically." Example Prt 6 Required setting dout Interval output mode

> dout Interval time ınŀ

#### Caution

The balance may not transmit the data completely at the specified interval times and baud rate. Set the baud rate higher.

# 9-6. Description Of The Item "Data Format"

#### A&D Standard Format

SIF HYPF N

This format is used when the peripheral equipment can receive the A&D format. If an AD-8121B printer is used, set the printer to MODE 1 or 2.

- This format consists of fifteen characters excluding the terminator.
- A header of two characters indicates the balance condition.
- The polarity sign is placed before the data with the leading zeros. If the data is zero, the plus sign is used.
- The unit, consisting of three characters, follows the data.
- Output example: S T, + 0 0 0 1. 2 7 8 \_ | c | t Unit **Terminator** Header Data Q T Stable header (Counting mode) T Stable header S Unstable header O L Overload header W T + | 1 Print example: С Header Data Unit Terminator


W T Header for modes other than the counting mode Q | T | Header for the counting mode

#### **DP (Dump Print) Format**

SIF EYPE I

This format is suitable for the peripheral equipment that prints the received data as is. When an AD-8121B printer is used, set the printer to MODE 3.


- This format consists of sixteen characters excluding the terminator.
- A header of two characters indicates the balance condition. No overload header is used.
- The polarity sign is placed before the data, with spaces in place of leading zeros, if the data is not zero or overloaded.
- The unit, consisting of three characters, follows the data.



KF Format 5 iF LYPE 2

This is the Karl-Fischer moisture meter format and is used when the peripheral equipment can only communicate using this format.


- □ This format consists of fourteen characters excluding the terminator.
- This format has no header characters.
- The polarity sign is placed before the data, with spaces in place of leading zeros, if the data is not zero or overloaded.
- This format outputs the unit only for a stable value.



MT Format 5 iF LYPE 3

This format is used when the peripheral equipment of other manufacturer is connected. Please note that the connection is not guaranteed.

- A header of two characters indicates the balance condition.
- The polarity sign is used only for negative data.
- The weight data uses spaces in place of the leading zeros.
- The character length of this format changes dependent upon the unit



#### **NU (Numerical) Format**

5 if LYPE 4

This format outputs only numerical data.

- □ This format consists of nine characters excluding the terminator.
- □ The polarity sign is placed before the data with the leading zeros. If the data is zero, the plus sign is used.

CSV Format 5 iF LYPE 5

- Separates the data of A&D standard format and the unit by a comma (,).
- Outputs the unit even when the data is overloaded.

$$\underbrace{\text{ST,+0001.278}}_{\text{Weight data}}, \ \_c \ \text{t}$$

S	Т	,	+	0	0	0	1		2	7	8	,		С	t	$C_R$	L <sub>F</sub>		
0	L	,	+	9	9	9	9	9	9	9	Е	+	1	9	,		С	t	C <sub>R</sub> L <sub>F</sub>

# 9-7. Data Format Examples

#### Stable

° 0.127 =t

A&D	S	Т	,	+	0	0	0	0		1	2	7		С	t	$c_R$	LF	
DP	W	Т		ш				+	0		1	2	7		С	t	$C_{R}$	LF
KF	+					0		1	2	7		С	t		$c_R$	LF		
MT		]	]	J	]	J	J	0		1	2	7	С	t	$c_R$	LF		
NU	+	0	0	0	0		1	2	7	$c_{R}$	ЧF							

#### **Unstable**

-18.369 <sup>= t</sup>

A&D	U	S	,	-	0	0	1	8		3	6	9		С	t	$C_{R}$	LF	
DP	U	ഗ	]	J	]	]	1	1	8	•	თ	6	9	]	O	t	$c_{R}$	$L_F$
KF	-	]	]	]	1	8		3	6	9					$c_R$	LF		
ΜT	S	О	]	]			1	8		3	6	9	С	t	$c_R$	LF		
NU	_	0	0	1	8		3	6	9	$C_{R}$	LF							

#### Overload Positive error

E st

A&D	0	L	,	+	9	9	9	9	9	9	9	Е	+	1	9	$c_R$	LF	
DP	J	I	J	J	]	]	J	J		J	J	J	J		J	J	$c_R$	LF
KF				]	J	]	Н	J							$c_R$	LЕ		
ΜT	S	I	+	$c_{R}$	ታ													
NU	+	9	9	9	9	9	9	9	9	$C_{R}$	ЧF							

#### Overload Negative error

-E =t

A&D	0	L	,	-	9	9	9	9	9	9	9	Е	+	1	9	$C_{R}$	LF	
DP				ш				-	Е	ш	ш			]		ш	$c_R$	LF
DP KF							L			ш					$C_{R}$	LF		
MT	S	I	-	$c_R$	LF													
NU	_	9	9	9	9	9	9	9	9	CR	LF							

∟ Space, ASCII 20h

c<sub>R</sub> Carriage Return, ASCII 0Dh

L<sub>F</sub> Line Feed, ASCII 0Ah

Units	Display	A&D	D.P.	KF	MT
g	g	g		g	ഥ g
Counting mode	PC5	⊔ P C	⊔ P C	□ p c s	□ P C S
Precent mode	%	<u> </u>	<u> </u>		山 %
Ounce (Avoir)	<u> </u>	O Z	⊔ O Z	_ 0 Z _	□ O Z
Pound	L <b>b</b>	□Ib	⊔Ib	u I b u	⊔ I b
Pound Ounce	L OZ	⊔ O Z	_ O Z	_ 0 Z _	□ O Z
Troy Ounce	ūΖt	o z t	o z t	니 o z t	ப o z t
Metric Carat	Ľť	∟ c t	_ c t	∟ c t ∟	_ c t
Momme	Min	m o m	m o m		ഥ m o
Pennyweight	d## t	d w t	d w t	ᆸ d w t	ᆸᅦd w t
Grain	5N	□GN	□GN	□ g r □	∟ G N
Tael (HK general, Singapore)	EL	⊔ t I	山 t I	∟ t I s	山 t I
Tael (HK, jewelry)	ŁL	山 t I	山 t l	니 t l h	山 t l
Tael (Taiwan)	EL	山 t I	山 t l	∟ t l t	∟ t l
Tael (China)	EL	山 t I	山 t I	⊔ t I c	山 t l
Tola (India)	EoL	ப ப t	ഥ t	⊔ t O I	□ t
Messghal	MES	m e s	m e s	□MS□	山 m
Multi	MLT	MLT	MLT	□ M L T	<u></u> М L Т


<sup>∟</sup> Space, ASCII 20h

#### Note

□ When "Pound Ounce" is selected, the data is output with the unit of ounce (oz).

# 9-8. Weighing Units

With the FZ-CT /FX-CT series balance, the following weighing units and weighing modes are available:



A unit or mode can be selected and stored in the function table as described on page 13. If a weighing mode (or unit of weight) has been turned off, that mode or unit will be missing in the sequence. Tael has four varieties, one of which can be selected and installed at the factory. To select a unit or mode for weighing, press the MODE key.

For details about the units and modes, see the table below:

Name (unit, mode)	Display	Function table (Storing mode)	Conversion factor 1 g =		
Gram	g	g	1 g		
Counting mode	P[5	PC <b>5</b>			
Percent mode	%	%			
Ounce (Avoir)	% [][]		28.349523125 g		
Pound	Lb	Lb	453.59237 g		
Pound/Ounce	L 07	LŪ	1 lb = 16 oz 1 oz = 28.349523125 g		
Troy Ounce	ΩZt	ΩZt	31.1034768 g		
Metric Carat	Ľt	Ľί	0.2 g		
Momme	m <u>n</u> m	mшm	3.75 g		
Pennyweight	dwt	₫#4‡	1.55517384 g		
Grain (UK)	5N	5N	0.06479891 g		
Tael (HK general, Singapore)			37.7994 g		
Tael (HK jewelry)	<u>-</u>	L-1	37.429 g		
Tael (Taiwan)		EL	37.5 g		
Tael (China)			31.25 g		
Tola (India)	EoL	EoL	11.6638038 g		
Messghal	ME5	ME5	4.6875 g		
Programmable-unit (Multi-unit)	MLT	MLT			

The tables below indicate the weighing capacity and the minimum display for each unit, depending on the balance model.

		Capacity		
Unit	FZ-1200CT FX-1200CT	FZ-700CT FX-700CT	FZ-500CT FX-500CT	Minimum Display
Gram	252	152	102	0.001
Ounce (Avoir)	8.89	5.36	3.59	0.00005
Pound	0.555	0.335	0.224	0.000005
Pound/Ounce	0L 8.88 oz	0L 5.36 oz	0L 3.60 oz	0L 0.01 oz
Troy Ounce	8.10	4.88	3.27	0.00005
Metric Carat	1260	760	510	0.001
Momme	67.2	40.5	27.2	0.001
Pennyweight	162.0	97.7	65.5	0.001
Grain (UK)	3888	2345	1574	0.02
Tael (HK general, Singapore)	6.66	4.02	2.69	0.00005
Tael (HK jewelry)	6.73	4.06	2.72	0.00005
Tael (Taiwan)	6.72	4.05	2.72	0.00005
Tael (China)	8.06	4.86	3.26	0.00005
Tola (India)	21.6	13.0	8.74	0.0001
Messghal	53.7	32.4	21.7	0.0005

# 9-9. Storing Units

The units or modes can be selected and stored in the function table. The sequence of displaying the units or modes can be arranged in the function table so as to fit the frequency of use.

Select a unit or mode and arrange the sequence of display as follows:


- 1 Press and hold the SAMPLE key until b85Fnc of the function table is displayed.
- 2 Press the SAMPLE key several times to display Unit.
- 3 Press the PRINT key to enter the unit selection mode.
- 4 Specify a unit or mode in the order to be displayed using the following keys.

SAMPLE key ......To sequentially display the units.

RE-ZERO key......To specify a unit or mode.

The stabilization indicator appears when the displayed unit or mode is specified.

- 5 Press the PRINT key to store the units or modes. The balance displays
  - $\boxed{\mathcal{E}_{nd}}$  and then displays the next menu item of the function table.
- 6 Press the CAL key to exit the function table. Then the balance returns to the weighing mode with the unit specified first in step 4.



### **Notes**

- □ When the power is turned on, the scale displays the unit specified first in step 4.
- □ In the weighing mode, to select a unit or mode for weighing, press the MODE key.

# 9–10. Counting Mode (PCS)

This is the mode to determine the number of objects in a sample based on the standard sample unit mass. The unit mass means the mass of one piece of the sample. The smaller the variables in each sample unit mass are, the more accurate the counting will be. The balance is equipped with the Automatic Counting Accuracy Improvement (ACAI) function to improve the counting accuracy.

#### **Notes**

- For counting, use samples with a unit mass of at least ten or more times the minimum weighing value of the balance.
- □ If the sample unit mass variable is too large, it may cause a counting error.
- To improve the counting performance, use the ACAI function frequently or divide the samples into several groups and count each group.

# **Selecting The Counting Mode**

1 Press the MODE key to select F55 (counting mode).

## **Storing A Sample Unit Mass**

- 2 Press the SAMPLE key to enter the sample unit mass storing mode.
  - Even in the storing mode, pressing the MODE key will switch to the next mode.
- 3 To select the number of samples, press the SAMPLE key several times. It may be set to 10, 25, 50 or 100.


### Note

- A greater number of samples will yield more accurate counting result.
- Place a container on the weighing pan, if necessary.

  Press the RE-ZERO key to cancel the weight (tare).

  The number specified in step 3 appears.

  e.g.: 25 0 is displayed if 25 is selected in step 3.
- 5 Place the number of samples specified on the pan. In this example, 25 pieces.
- Wait for the stabilization indicator to turn on. Press the PRINT key to calculate and store the unit mass. The balance displays [-5-5-3] (counting mode) and is set to count samples with this unit mass. To improve the accuracy of the unit mass, go to step 8.



### The counting result

#### Notes

- If the balance judges that the mass of the samples is too light to aquire accurate weighing, it displays an error requiring the addition of more samples to the specified number. For example, 50 53 appears to require 25 more samples. Add 25 samples and press the PRINT key. When the unit mass is stored correctly, the balance proceeds to the counting mode.

- For a more accurate counting operation, use samples with a unit mass of at least ten or more times the minimum weighing value of the balance. For example, when a balance with a minimum weighing value of 0.0001 g is used, the recommended sample unit mass is 0.001 g or more.
- The sample unit mass stored, even if the AC adapter is removed, is maintained in non-volatile memory.

### **Counting Operation**


7 Place the samples to be counted on the pan.
While the stabilization indicator is on, pressing the PRINT key will output the weight value (number of objects), using the RS-232C serial interface.

#### Note

 Peripheral equipment, that is sold separately, such as a printer or a personal computer is required.

### **Counting Mode Using The ACAI Function**

The ACAI is a function that improves the accuracy of the unit mass automatically by increasing the number of samples and averaging the unit mass variable to minimize the weighing error, as the counting process proceeds.



- 8 If a few more samples are added, the processing indicator turns on. To prevent an error, add three or more. The processing indicator does not turn on if overloaded. Try to add the same number of samples as displayed.
- 9 The balance re-calculates the unit mass while the processing indicator is blinking. Do not touch the balance or samples on the pan until the processing indicator turns off.
- 10 Counting accuracy is improved when the processing indicator turns off.
  Each time the above operation is performed, a more accurate unit mass will be obtained. There is no definite upper limit of ACAI range for the number of samples exceeding 100. Try to add the same number of samples as displayed.
- 11 Remove all the samples used in ACAI and proceed with the counting operation using the improved unit mass.

#### Percent Mode (%) 9–11.

This is the mode to display the weight value in percentage compared to a 100% reference mass and is used for target weighing or checking the sample variable.

## Selecting The Percent Mode

Press the MODE key to select % (percent mode).

### Storing The 100% Reference Mass

- 2 Press the SAMPLE key to enter the 100% reference mass storing mode.
  - Even in the storing mode, pressing the MODE key will switch to the next mode.
- Place a container on the weighing pan, if necessary. Press the RE-ZERO key to cancel the weight (tare). The balance displays 100 0 %.
- Place the sample to be set as the 100% reference mass on the pan or in the container.
- Press the PRINT key to store the reference mass. The balance displays | 100,000 % |. (The decimal point position depends on the reference value.)

#### **Notes**


- If the balance judges that the mass of the sample is too light to be used as a reference, it displays La. Do not use the sample.
- The reference mass stored, even if the AC adapter is removed, is maintained in non-volatile memory
- 6 Remove the sample.


# Reading The Percentage

Place a sample to be compared to the reference mass on the pan. The displayed percentage is based on the 100% reference mass. While the stabilization indicator is on, pressing the PRINT key will output the weight value, using the RS-232C serial interface.

#### Note

Peripheral equipment, that is sold separately, such as a printer or a personal computer is required.





Weighing result compared to 100% reference mass

#### **Description Of The Item "Application Function"** 9–12.

# Capacity Indicator ( RPF | 1)

In the weighing mode, the indicator displays the weight data relative to the weighing capacity in percentage. (Zero = 0%, maximum capacity = 100%)

# Statistical Calculation Mode ( RPF 2)

The mode statistically calculates the weight data, and displays or outputs the results. For details, refer to "6-5 Statistical Calculation Mode".

# 9-13. Statistical Calculation Mode

The statistical calculation mode statistically calculates the weight data, and displays or outputs the results. To use the statistical calculation mode, set the "Application function ( $\Pi P = F \cap E$ )" in the function table to "E", as described below.

Statistical items available are number of data, sum, maximum, minimum, range (maximum-minimum), average, standard deviation, coefficient of variation and relative error. What statistical items to output can be selected from the four modes in the function table.

- The wrong data input can be canceled by the key operation, if immediately after the input.
- Turning the balance off will delete the statistical data.
- The standard deviation, coefficient of variation and relative error are obtained by the equation below:

Standard deviation = 
$$\sqrt{\frac{\sum (X_i)^2 - (\sum X_i)^2}{N \cdot (N-1)}}$$
 where Xi is the i-th weight data, N is number of data.

Coefficient of variation (CV)=  $\frac{\text{Standard deviation}}{\text{Average}} \times 100 \, (\%)$ 

Relative error of maximum value =  $\frac{\text{Maximum value - Average}}{\text{Average}} \times 100 \, (\%)$ 

Relative error of minimum value =  $\frac{\text{Maximum value - Average}}{\text{Average}} \times 100 \, (\%)$ 

#### Note

 When data with the minimum weighing value turned off are used in the statistical calculation, results will be displayed with the minimum weighing value turned off. (The minimum weighing value will be rounded off.)

### 9.13.1. Getting Started


# Switching To The Statistical Function Mode (Changing The Function Table)

Press and hold the SAMPLE key until 68SEnc 1/10d Press and hold SAMPLE the function table is displayed, then release the key. LASFne Press the SAMPLE key several times to display 2 1/10d Press several times AP Foc SAMPLE AP Foc Press the PRINT key to display ⊙ 4 Press the RE-ZERO key several times to display Seft To select statistical items to output, go to step 5. Mari To store the statistical function mode setting, go to step 7. To disable the statistical calculation mode, press the →0← RE-ZERO Press several times RE-ZERO key to select RPF Na--APF 5ERtTo store the setting To select statistical go to step 7. go to step5.

# Selecting The Statistical Items To Output

- 5 Press the SAMPLE key to display رقيق المجاهة (SAMPLE) key to display
- 6 Press the RE-ZERO key to select the output items.
  In the example, Star is selected to output the number of data, sum, maximum, minimum, range (maximum-minimum) and average.

Parameter	Description
<b>■</b> □	Number of data, Sum
1	Number of data, Sum, Maximum, Minimum,
,	Range (maximum – minimum), average
	Number of data, Sum, Maximum, Minimum,
2	Range (maximum – minimum), Average, Standard deviation,
	Coefficient of variation
	Number of data, Sum, Maximum, Minimum,
] =	Range (maximum – minimum), Average, Standard deviation,
	Coefficient of variation, Relative error of maximum value,
	Relative error of minimum value




7 Press the PRINT key to store the setting.

#### Note

- The statistical calculation mode will be enabled using the ON:OFF key or when the power is turned on.
- 8 Press the CAL key to return to the weighing mode.

## **Selecting The Unit**

Press the  $\boxed{\text{MODE}}$  key to select the unit to be used for the statistical calculation mode. In the example shown at the right, carat ( $\not = t$  ) is selected.



#### Notes

- Selecting the unit using the MODE key is not available after the data is entered. In this case, clear the data as described on page 43 and select the unit using the MODE key.
- □ When the unit used for the statistical calculation mode is to be enabled upon power-on, select the unit in "Unit ( $Un \ i E$ )" of the function table beforehand.

# 9.13.2. Using The Statistical Calculation Mode

# **Entering Data For Statistical Calculation**

Use the following keys to operate the statistical calculation mode.

- MODE key ...... When the data is entered, moves between the displaying items (weighing mode, statistical results and data operation) each time the key is pressed.
  - When no data has been entered, selects the unit.

SAMPLE key..... Turns the minimum weighing value on or off, in the weighing mode.


RE-ZERO key  $\cdots$  Sets the display to zero in the weighing mode.

PRINT key ...... Outputs the data number and the weight data and includes the weight data to statistical calculation in the weighing mode. (Output is not in the data format specified in the function table because of the data number added.)

 Outputs the statistical results while the statistical results are displayed. (Output is not in the data format specified in the function table.)

CAL key..... Returns to the weighing mode.

- 1 Press the RE-ZERO key to set the display to zero.
- 2 Place the sample on the weighing pan and wait for the stabilization indicator to turn on.
- 3 Press the PRINT key to add the data displayed to statistical calculation. The number of data on the upper left of the display increases by 1.
- 4 Repeat steps 1 to 3 for each weighing.




## **Outputting The Statistical Results**

- 5 Each time the MODE key is pressed, the display changes: the results as selected in "Statistical function mode output items (5ŁRF)", ELERR and ERNEEL.
- 6 While the results are displayed, press the PRINT key to output the results.

#### Notes

- □ When the number of data is 1, the coefficient of variation is displayed as -----.
- □ When the average is 0, the coefficient of variation is displayed as -----.
- Statistical items are indicated on the upper left of the display using the following symbols.

Symbol	Statistical item
Sun	Sum
āRH	Maximum
_ _ iu	Minimum
г	Range (Maximum – minimum)
RūE	Average
5d	Standard deviation
Εū	Coefficient of variation
āRH%	Relative error of maximum value
- n In%	Relative error of minimum value



# Output example Function table parameter


N SUM	10	
	100.0000g	
MAX	10.5000g	]   2
MIN	9.5000g	3
R	1.0000g	
AVE	10.0000g	
SD	0.280 g	
CV	2.80 %	
MAX%	5.00 %	
MIN%	5.00 %	<b>]</b>

# 9.13.3. Deleting The Latest Data

When the wrong data is entered, it can be deleted and excluded from statistical calculation. Only the latest data can be deleted.

- 1 In the weighing mode, press the MODE key several times to display [FRIEFL].
- 2 Press the PRINT key to display [[RN no.]
- 3 Press the RE-ZERO key to display [ AN Lo
- 4 Press the PRINT key to delete the latest data and exclude it from statistical calculation.


The number of data decreases by 1 when the balance returns to the weighing mode.




# 9.13.4. Clearing The Statistical Data

All the statistical data will be deleted and the number of data will be 0 (zero).

- 1 In the weighing mode, press the MODE key several times to display [1.588].
- 2 Press the PRINT key to display [[Lr na].
- 3 Press the RE-ZERO key to display [[Lr [[]]].
- 4 Press the PRINT key to delete the statistical data. The number of data becomes 0 (zero) when the balance returns to the weighing mode.





# 9–15. Clock And Calendar Function (Only For The FZ-CT Series)

The **FZ-CT** series balance is equipped with a clock and calendar function. When the "GLP output  $(\sqrt{n} + a)$ " parameter is set to "/" or "/" and the "Time/Date output (5 - b/a)" parameter is set to "/", "/" or "/", the time and date are added to the output data. Set or confirm the time and date as follows:

### **Operation**

- 1 Press and hold the SAMPLE key until b85Fnc of the function table is displayed.
- 2 Press the SAMPLE key several times to display [[L Rd]].
- 3 Press the PRINT key. The balance enters the mode to confirm or set the time and date.

# Confirming the date

- 4 The current date is displayed with all the digits blinking.
  - When the date is correct and the operation is to be finished, press the CAL key and go to step 8.
  - When the date is correct and the time is to be confirmed, press the SAMPLE key and go to step 6.
  - □ When the date is not correct and is to be changed, press the RE-ZERO key and go to step 5.

#### Note

The year is expressed using a two-digit format. For example, the year 2010 is expressed as "10".

# Setting the date (with part of the digits blinking)

5 Set the date using the following keys.

SAMPLE key To select the digits to change the value. The

selected digits blink.

RE-ZERO key To increase the value by one.

MODE key To decrease the value by one.


PRINT key To store the new setting, display End

and go to step 6.

CAL key To cancel the new setting and go to step 6.

# Confirming the time

- 6 The current time is displayed with all the digits blinking.
  - When the time is correct and the date does not need to be confirmed, press the CAL key and go to step 8.
  - When the time is correct and the date is to be confirmed, press the SAMPLE key and go back to step 4.
  - When the time is not correct and is to be changed, press
     the RE-ZERO key and go to step 7.



### Setting the time (with part of the digits blinking)

7 Set the time in a 24-hour format using the following keys.

SAMPLE key To select the digits to change the value.

The selected digits blink.

RE-ZERO key To increase the value by one.

MODE key To decrease the value by one.

PRINT key To store the new setting, display End

and go to step 8.


CAL key To cancel the new setting and go to step 4.

## Quitting the operation

8 The balance displays the next menu of the function table. Press the CAL key to exit the clock and calendar function and return to the weighing mode.

### **Note**

 Do not enter invalid values such as a non-existing date when setting the time and date.



# 10.ID Number And GLP Report

- The ID number is used to identify the balance when Good Laboratory Practice (GLP) or Good Manufacturing Practice (GMP) is used.
- □ The ID number is maintained in non-volatile memory even if the AC adapter is removed.
- The output format for GLP/GMP compliant report is selected at "GLP output ( $nF_0$ )" of the function table and can be output to a personal computer or printer using the RS-232C serial interface.
- The GLP/GMP compliant report includes the balance manufacturer, model, serial number, ID number and space for signature for weight data, and the weight used and results for calibration or calibration test data.
  - When the AD-8121B printer is used, the date and time can be printed using the AD-8121B clock and calendar function. In this case, set the "GLP output ( mFa)" parameter to " l".
- □ The balance can output the following:
  - "Calibration report" of the calibration, using the internal mass (Calibration due to changes in temperature and one-touch calibration.)
  - "Calibration report" of the calibration, using an external weight.
  - "Calibration test report" of the calibration test, using an external weight and using the internal mass.
  - "Title block" and "End block" for the weight data.

# 10–1. Setting The ID Number

- 1 Press and hold the SAMPLE key until bRSFnc of the function table is displayed.
- 2 Press the SAMPLE key several times to display d.
- 3 Press the PRINT key. Set the ID number using the following keys.

SAMPLE key To select the digit to change the value. The selected digit blinks.

To set the character of the digit selected. Refer to the display character set shown below.

 PRINT
 key
 To store the new ID number and display
 RP Fnc

 CAL
 key
 To cancel the new ID number and display
 RP Fnc

4 Press the CAL key to return to the weighing mode.

### Display character set

RE-ZERO key, MODE key

0	1	2	3	4	5	6	7	8	9		Г	Α	В	С	D	Ε	F	G	Н	I	J	K	L	М	N	0	Р	Q	R	s	Т	U	V	W	X	Υ	z
	11		3	닉	5	6	7	8	9	1		Я	Ь	1	Ы	E	Ë	ū	Н	1	П	K	L	M	N	O	₽		R	בה	Ŀ	Ш	1/	H	> <	닠	7
<			-[N	10[	DΕ	key						Spa	асе																F	RE-	ZE	RO	ke	у -			$\overline{\rightarrow}$

# 10–2. **GLP Report**

Set the following parameters to output the GLP/GMP compliant report.

- To print the report, set the "GLP output ( InFa )" parameter to " /", set the "Data output pause ( PUSE )" parameter to " /" and use MODE 3 of the AD-8121B printer. For details on using the printer, refer to "16-1 Connection to the AD-8121B Printer".
- □ To output the report to a personal computer using the RS-232C serial interface, set the "GLP output ( ¬¬F¬□ )" parameter to "¬□".

#### Note


□ If the time and date are not correct, adjust the AD-8121B clock and calendar (FZ-CT series) or use "[L Rdd" of the function table to adjust the time and date (FZ-CT series).

### Calibration Report Using An Internal Mass (Only For The FZ-CT Series)

When the setting is " InFa !":

When the setting is " ¬¬F¬ ¬?":

### AD-8121 format



## Calibration Test Report Using The Internal Mass (Only For The FZ-CT Series)

(Calibration test does not perform calibration.)


LF Line feed, ASCII 0Ah

LF Line feed, ASCII 0Ah


When the setting is " InFa !":

When the setting is " ¬¬F¬¬ ¬?":

### AD-8121 format



### Calibration Report Using An External Weight



With FZ-CT series, the time and date are added to output data.

# Calibration Test Report Using An External Weight

(Calibration test does not perform calibration.)

LF Line feed, ASCII 0Ah


When the setting is " nFa | l":

LF Line feed, ASCII 0Ah

AD-8121 format

When the setting is " ¬¬F¬¬¬ ¬?":

General data format



\* With FZ-CT series, the time and date are added to output data.

### Title Block And End Block


When a weight value is recorded as the GLP data, "Title block" and "End block" are inserted at the beginning and at the end of a group of weight values, in the GLP report.

#### Note

□ To output the report to an AD-8121B, use MODE 3 of the AD-8121B.

### Operation

- 1 With the weight data displayed, press and hold the PRINT key until 5 kg unt
- 2 The weight data is output according to the parameter setting of the data output mode.
- 3 Press and hold the PRINT key until ΓΕςΕπά is displayed. The "End block" is output.




# 11.Underhook

The underhook can be used for magnetic materials or density measurement. The built-in underhook is revealed by removing the plastic cap on the bottom of the balance.

- Underhook

Cap

Use the underhook as shown below.



### Caution

Do not apply excessive force to the underhook.

When not in use, attach the plastic cap to prevent dust from getting into the balance.

Do not push the underhook upward.

# 12.Programmable-Unit

This is a programmable unit conversion function. It multiplies the weight data in grams by an arbitrary coefficient set in the function table and displays the result.

The coefficient must be within the range between the minimum and maximum shown below. If the coefficient set is beyond the range, an error is displayed and the balance returns to the coefficient setting mode, prompting to enter an appropriate value. A coefficient of 1 was set at the factory.

	Model		Minimum coefficient	Maximum coefficient
FZ-1200CT	FZ-700CT	FZ-500CT	0.000001	1000
FX-1200CT	FX-700CT	FX-500CT	0.00001	1000

## Operation

- 1 Press and hold the SAMPLE key until bR5Fnc of the function table is displayed.
- 2 Press the SAMPLE key several times to display MLT.
- 3 Press the PRINT key. The balance enters the mode to confirm or set the coefficient.

# Confirming the coefficient

- 4 The current coefficient is displayed with the first digit blinking.
  - □ When it is not to be changed, press the CAL key and go to step 6.
  - When it is to be changed, press the RE-ZERO key and go to step 5.

# Setting the coefficient

5 Set the coefficient using the following keys.

SAMPLE key ......To select a digit to change the value.

The selected digit blinks.

RE-ZERO key......To change the value.

MODE key .....To change the decimal point position.

Each time the switch is pressed, the decimal

point position changes as follows:

PRINT key......To store the new setting, display End and go to step 6.

CAL key.....To cancel the new setting and go to step 6.

# Quitting the operation

6 The balance displays Lin L. Press the CAL key to exit the programmable-unit function and return to the weighing mode.

# 

# Using the function

Press the MODE key to select the programmable-unit (no display on the unit section). Perform weighing as described in "5–1. Basic Operation (Gram Mode)". After weighing, the balance displays the result (weight data in grams x coefficient).

# 13. Password Function

The password function can be specified for an authorized user by using a password stored in the balance by manager. While the password function is used, the balance can not be weighed without inputting the correct password.

Password: Four digits (4 x 4 x 4 x 4 = 256 outcomes) using MODE, SAMPLE, PRINT and RE-ZERO keys.

At factory setting, the password function is not used.


The management of the password function is performed in "Password function" of "9–3. Details Of The Function Table".

# 13-1. Using Password Function

The password function can be switched between "Use" or "Not used" at "Password function ( $L \square \varepsilon F \cap \varepsilon$ )" of "9–3. Details Of The Function Table".

- 1 Press and hold the SAMPLE key in the weighing mode.

  BRSFnc is displayed.
- 2 Press the SAMPLE key several times until LacFac is displayed.
- 3 Press the PRINT key to display PASS To cancel the operation, press the CAL key.
- 4 Press the RE-ZERO key to display PR55 XX
- 5 Press the PRINT key to display Sur E No" where "No" is selected.
- 6 "YES" and "No" can be switched with the RE-ZERO key. Press the RE-ZERO key to display 5ur E No W where "YES" is selected.
- 7 Press the PRINT key to store the new setting. The password function can now be used.
  When turning on the balance, the password input is displayed.



# 13-2. Changing Password


The password can be changed at "Password (Lac na)" of "9-3. Details Of The Function Table".

- 1 Press and hold the SAMPLE key in the weighing mode.

  \$\begin{align\*} \begin{align\*} \begin{ali
- 2 Press the SAMPLE key several times until Lac na is displayed.
- 3 Press the PRINT key to display the current password.
  At factory settings, the password is 7777. (the RE-ZERO key, four times)
- 4 Set the new password using the following keys. The balance will turn automatically after no operation for ten minutes.

MODE key Character	M
SAMPLE key Character	5
PRINT key Character	P
RE-ZERO key Character	7
CAL key Back key	

- 5 Input four characters of the new password using these keys.
- 6 The balance displays 5ur E No where "No" is blinking when "No" is selected.
- 7 Press the RE-ZERO key to display 5ur E No M where "YES" is blinking when "YES" is selected.



### **Note**

If the correct password is missing, the balance can not be used.
 Record the correct password after changing.

# 13-3. Inputting Password When Turning On The Balance

When turning on the balance using the password function, the password input is displayed before weighing.

1	Turn on the balance using the ON:OFF key.	-	
2	$PR55$ is displayed and ${}_{\#}$ of the password input is displayed.	ON:OFF	
3	Input the four characters of the correct password using the following keys.  MODE key	#5\\\\	After no operation for 10 minutes
4	When the password is correct, the balance displays weighing value after full-segment display.  When the password is incorrect, FR L is displayed and the buzzer sounds three times.	Correct  Correct  Weighing display	F F I I L

# 13-4. Missing Password

If the correct password is missing, the balance can not be used. Contact your local A&D dealer to reset the password to factory settings.

# 14.RS-232C Serial Interface

The balance is a Data Communication Equipment (DCE) device. Connect the balance to a personal computer (DTE) using a straight through cable.

Transmission system : EIA RS-232C (D-Sub 9-pin, female connector)

Transmission form : Asynchronous, bi-directional, half duplex

Transmission rate : Approx. 10 or 5 times/second (same as data refresh rate)
Data format : Baud rate : 600, 1200, 2400, 4800, 9600, 19200 bps

Data bits : 7 or 8 bits


Parity : Even, Odd (Data bits 7 bits)


None (Data bits 8 bits)

D-Sub 9-Pin Assignments

Pin No.	Signal name	Direction	Description
1	-	-	No connection
2	TXD	Output	Transmit data
3	RXD	Input	Receive data
4	-	-	No connection
5	SG	-	Signal ground
6	DSR	Output	Data set ready
7	RTS	Input	Request to send
8	CTS	Output	Clear to send
9	-	-	No connection

Signal names of the balance side are the same as the DTE side with TXD and RXD reversed.





# 15. Connection To Peripheral Equipment

# 15–1. Connection To The AD–8121B Printer

Set the following parameters to use the AD-8121B printer.

Example of use	AD-8121B mode setting
To print A&D standard format weight data, using the FZ-CT /FX-CT PRINT key or auto print mode. (The time and date can be added.)	MODE 1
To print A&D standard format weight data, using the AD-8121B DATA key or AD-8121B built-in timer. (The time and date can be added.) To print, using the AD-8121B chart printing function.	MODE 2
To print the FZ-CT / FX-CT statistical data.	MODE 3
To print GLP output.	MODE 3

Class	Item and Parameter	Factory setting	AD-8121B MODE 1	AD-8121B MODE 2	AD-8121B MODE 3
dout	Pr Ł Data output mode	0	0,1,2,4,5 *1	3	0,1,2,4,5 *1
Data output	PUSE Data output pause	0	0	0	0,1 *2
	6년5 Baud rate	2	2	2	2
5 ,F Serial	したりに Data bit, parity bit	0	0	0	0
interface	ErLF Terminator	0	0	0	0
	는	0	0	0	1

<sup>\*1</sup> Set appropriate parameters for "PP-P (Auto print polarity)" and "PP-b (Auto print difference) when auto print mode A or B ( $Prb \mid r$  or  $Prb \mid r$ ) is selected. Set the AD-8121B DIP switch No.3 to ON when unstable data is printed with " $Prb \mid r$ ".

### **Notes**

Refer to "10-2. GLP Report" for print samples. Settings of AD-8121B DIP switches

MODE	AD-8121B DIP switch	Description
MODE 1	D1234	Print at receiving data.
IVIODE I		Standard mode, statistical calculation mode
MODE 2	O1234	Print by the AD-8121B DATA key operation or AD-8121B built-in timer.
MODE 2		Standard mode, interval mode, chart mode
MODE 3	D1 2 3 4	Print at receiving data.
IVIODE 3		Dump print mode

DIP switch No.3: Handling unstable data

ON Print

OFF Not printed

D 1 2 3 4

Set the DIP switch No.4 to OFF.






<sup>\*2</sup> Set / when multiple lines are printed in the **FZ-CT** / **FX-CT** statistical calculation mode.

# 15–2. Connection To A Computer

The balance can be connected to a personal computer using the RS-232C serial interface. As an option, the FXi-02 USB interface is available to transmit the balance data to a personal computer.

The balance is a DCE device. Use a straight through cable when connecting to a personal computer. A commercially available modem cable can be used. When purchasing, confirm the following specifications.



# 15-3. Using Windows Communication Tools (WinCT)

When Windows is used as an operating system in a personal computer, the WinCT software, that is downloaded from the A&D website, can be used to transmit the weight data to the personal computer. The WinCT has three communication methods: "RsCom", "RsKey" and "RsWeight".

### **Rscom**

- Can transmit commands to control the balance.
- Can make bi-directional communication between the balance and a personal computer using the RS-232C serial interface.
- Can display or store the data using a text file format. Can also print the data using a printer connected to a personal computer.
- When several ports of a personal computer have balances connected, can communicate with each balance simultaneously.
- Can share a personal computer with other application software.

### **Rskey**

- Can transmit the weight data output from the balance directly to other application software such as Microsoft Excel.
- Can be used with most application software.

## Rsweight

- Can retrieve the weight data from the balance and display the data as a graph in real time.
- Can calculate and display the maximum, minimum, average, standard deviation and coefficient of variation.

# Using The WinCT Software, The Balance Can Do The Following:

- 1 Analyzing the weight data and the statistics with "RsKey"
  - The weight data can be input directly into an Excel worksheet. Then, Excel can analyze the data to obtain sum, average, standard deviation, maximum and minimum value, and display them as a graph.
- 2 Controlling the balance using commands from a personal computer By using "RsCom", the personal computer sends commands such as "re-zero" or "send weight data" to the balance and controls the balance.
- 3 Printing the balance GLP report using your printer
  The balance GLP report can be printed using a printer connected to a personal computer.
- 4 Receiving weight data at a certain interval The weight data can be received at a certain interval and data characteristic with elapsed time can be obtained.

Using a personal computer as an external indicator
With the "RsKey" test mode function, a personal computer can be used as an external weight indicator for the balance. (To do this, set the balance data output mode to stream mode.)

# 16.Commands

# 16-1. Command List

### Note

A command has a terminator added, that is specified using the "Terminator ( $[\Gamma LF]$ )" parameter of "Serial interface ( $[\Gamma LF]$ )" in the function table, and is sent to the balance.

Commands To Query Weight Data

	, - 5
С	Cancels the S or SIR command.
Q	Requests the weight data immediately.
S	Requests the weight data when stabilized.
SI	Requests the weight data immediately.
SIR	Requests the weight data continuously.
ESCP	Requests the weight data when stabilized.

Note: The "Q" and "SI" commands, the "S" and "ESCP" commands behave the same.

### **Commands To Control The Balance**

	of the balance
CAL	Same as the CAL key.
EXC *	Calibration using an external weight
OFF	Turns the display off.
ON	Turns the display on.
Р	Same as the ON:OFF key
PRT	Same as the PRINT key
R	Same as the RE-ZERO key
SMP	Same as the SAMPLE key.
T	Tare key
Z	Same as the RE-ZERO key
ESCT	Same as the RE-ZERO key
U	Same as the MODE key
?ID	Requests the ID number.
?SN	Requests the serial number.
?TN	Requests the model name.
?PT	Requests the tare weight.
PT: ***.***_ct	Changes the tare weight.
	The unit added is the current weighing unit in A&D standard format.

### Note:

□ The "R", "Z" and "ESCT" commands behave the same.

 $^{\text{E}}\text{S}_{\mathbb{C}}$ : ASCII code 1Bh

\* : Only for the FZ-CT series

# 16-2. Acknowledge Code And Error Codes

When the "AK, Error code ( $\mathcal{E}_{\Gamma}\mathcal{E}_{d}$ )" parameter of "Serial interface ( $\mathcal{E}_{\Gamma}\mathcal{E}_{d}$ )" is set to " $\mathcal{E}_{\Gamma}\mathcal{E}_{d}$ )" the balance outputs <AK> code or an error code for each command as follows:

- <AK> (06h) Acknowledge in ASCII code.
- When the balance receives a command to request data and can not process it, the balance transmits an error code (EC, Exx).
  - When the balance receives a command to request data and can process it, the balance outputs the data.
- When the balance receives a command to control the balance and can not process it, the balance transmits an error code (EC, Exx).
  - When the balance receives a command to control the balance and can process it, the balance transmits the acknowledge code.


Among commands to control the balance, the following transmit the acknowledge code both when the balance receives the command and when the balance has accomplished the command. If the command can not be processed properly, the balance transmits an error code (EC, Exx). This error can be released using the CAL command.

CAL command (When performing the calibration using a calibration weight)


ON command, P command

R command, Z command, T command (When setting the display to zero)


## R command example



## T command example



## PT command example



□ When a communication error has occurred due to external noise, or a parity error has occurred due to transmission error, the balance transmits an error code. In this case, send the command again.

# 16-3. Settings Related To RS-232C

Concerning the RS-232C, the balance has two functions: "Data output (daub)" and "Serial interface ( $5 \, F$ )". Set each function as necessary.

# 17. Maintenance

 Do not use organic solvents to clean the balance. Clean the balance with a lint free cloth that is moistened with warm water and a mild detergent.

Balance main unit.	Use a soft, lint free cloth that is moistened with a mild detergent to clean.
Weighing pan	The edge of the weighing pan is sharp. Use caution when cleaning the pan.
Breeze break (standard accessory)	An antistatic treatment has been applied to the breeze break components. Use a soft, dry, lint free cloth to clean. Cleaning with a cloth that is moistened with water or a mild detergent repetitively, or washing with water, may degrade the antistatic effect.

- □ Do not disassemble the balance. Contact the local A&D dealer if the balance needs service or repair.
- Use the original packing material for transportation.

# 18. Troubleshooting

# 18-1. Checking The Balance Performance And Environment

The balance is a precision instrument. When the operating environment or the operating method is inadequate, correct weighing can not be performed. Place a sample on the pan and remove it, and repeat this several times. If the balance seems to have a problem with repeatability or to perform improperly, check as described below. Also, visit our website, http://www.aadd.jp/, for "Information Library" and "FAQ". If improper performance persists after checking, contact the local A&D dealer for repair.

## Checking That The Balance Performs Properly

- Check the balance repeatability using an external weight. Be sure to place the weight in the center of the weighing pan.
- Check the balance repeatability, linearity and calibrated value using external weights with a known value.

# Checking That The Operating Environment Or Weighing Method Is Proper Operating environment

- Is the weighing table solid enough?
- □ Is the balance level? Refer to "3-1. Before Use".
- Is the operating environment free from vibration and drafts? Has the breeze break been installed?
- Is there a strong electrical or magnetic noise source such as a motor near the balance?

### Weighing method

- Is the weighing pan installed correctly?
- Is the RE-ZERO key pressed before placing a sample on the weighing pan?
- Is the sample placed in the center of the weighing pan?
- Has the balance been calibrated using an external weight? (Or one-touch calibration, only for the FZ-CT series)
- Has the balance been warmed up for 30 minutes before weighing?

## Sample and container

- Has the sample absorbed or lost moisture due to the ambient conditions such as temperature and humidity?
- Has the temperature of the container been allowed to equalize to the ambient temperature? Refer to "3-2. During Use".
- □ Is the sample charged with static electricity? Refer to "3–2. During Use".
  - The FZ-CT / FX-CT series balances are prone to be charged with static electricity when the relative humidity is low.
- □ Is the sample of magnetic material such as iron? Caution is required for weighing magnetic materials. Refer to "3–2. During Use".

18-2. <b>E</b>	rror Codes	
Display	Error code	Description
		Stability error The balance can not stabilize due to an environmental problem. Prevent
Error !		vibration, drafts, temperature changes, static electricity and magnetic
	EC, E11	fields. Refer to "3. Precautions" for details on the operating environment
		and "6. Response Adjustment" about adapting the balance to the
		environment. To return to the weighing mode, press the CAL key.
Error2		Out of range error
		The value entered is beyond the settable range. Re-enter the value.
,		Internal mass error (Only for the FZ-CT series)
Error7	EC, E17	The internal mass application mechanism does not function properly.
		Perform the weighing operation from the beginning again.
		Calibration weight error
EAL E		The calibration weight is too heavy. Confirm that the weighing pan is
	EC, E20	properly installed. Confirm the calibration weight value.
		Press the CAL key to return to the weighing mode.
		Calibration weight error
-EAL E		The calibration weight is too light. Confirm that the weighing pan is
	EC, E21	properly installed. Confirm the calibration weight value.
		Press the CAL key to return to the weighing mode.
		Overload error
E		A sample beyond the balance weighing capacity has been placed on
		the pan. Remove the sample from the pan.
		Weighing pan error
		The weight value is too light. Confirm that the weighing pan and pan
- E		support are properly installed. Press the ON:OFF key two times to
		return to the weighing mode. If the error still persists, calibrate the
		balance.
		Sample mass error
Lo		The balance can not store the sample for the counting mode or for the
		percent mode because it is too light. Use a sample that is heavier.
<b>-</b> , -		Unit mass error
25 - *	<b>[</b> 5]	The sample unit mass for the counting mode is too light. Storing and
F 51 -		using it for counting will cause a counting error. Add samples to reach
50 - *	<b>I</b> 5	the specified number and press the PRINT key.
155		Pressing the PRINT key without adding samples will shift the balance
100 - *	<b>[5</b> ]	to the counting mode. But, to acquire accurate weighing, be sure to add
		samples.
Error0		Balance internal error
		If this error appears persistently, contact the local A&D dealer.
		Battery error for clock  The health is hotton in average for the clock was lost completely. After
		The back-up battery power for the clock was lost completely. After
rtc PF		pressing one of the keys, adjust the date and time. Even if the back-up
		battery power for the clock was lost and the balance is on or in the
		standby state, the balance performance will be normal. When this error
		appears frequently, repair is required.

Display	Error code	Description
	EC, E00	Communications error A protocol error occurred in communications. Confirm the format, baud
	·	rate and parity.
	EC, E01	Undefined command error
		An undefined command was received. Confirm the command.
	EC, E02	Not ready A received command can not be processed. e.g. The balance received a Q command, but not in the weighing mode.
		e.g. The balance received a Q command while processing a RE-ZERO command. Adjust the delay time to transmit a command.
	EC, E03	Timeout error  If the timeout parameter is set to "t-up", the balance did not receive the next character of a command within the time limit of one second.  Confirm the communication.
	EC, E04	Excess characters error The balance received excessive characters in a command. Confirm the command.
	EC, E06	Format error A command includes incorrect data. e.g. The data is numerically incorrect. Confirm the command.
	EC, E07	Parameter setting error The received data exceeds the range that the balance can accept. Confirm the parameter range of the command.
Other errors		If the errors described above can not be released or other errors are displayed, contact the local A&D dealer.

# 18-3. Other Display



The display ( PR55 " ) of the password input is displayed when using the password function and turning on the balance.

Enter a correct password. Refer to "13. Password Function".

# 18-4. Asking For Repair

If the balance needs service or repair, contact the local A&D dealer.

The balance is a precision instrument. Use much care when handling the balance and observe the following when transporting the balance.

- Use the original packing material.
- Remove the weighing pan from the main unit

# 19. Options And Peripheral Equipments

### **Note**

□ The FXi-02, FXi-08 and FXi-09 can not be used at the same time.

### FXi-02 USB interface

(Installed in the balance, Applicable OS: Windows 98 OSR2 or later)

- Used to transmit the balance weight data (numerical value only) uni-directionally to a personal computer via USB.
- Can transmit the balance weight data (numerical value only) directly to application software such as Microsoft Excel, Word and memo pad.
- Driver installation is not necessary.

### **Note**

 To perform bi-directional communication using WinCT or input statistical data or GLP data to a personal computer, use the USB converter (AX-USB-9P-EX)

FXi-08 Example of use



### FX*i*-08 Ethernet Interface

- Used to connect the balance to a LAN.
- The "WinCT-Plus" data communication software is provided as a standard accessory and can perform the following.
  - Acquire data from multiple balances connected to a LAN.
  - Control these balances with commands.
  - Acquire data transmitted from balances.
     Example: When pressing the PRINT key of the balance, data is output and is acquired by the computer.
  - The stored data can be used with Microsoft Excel (if installed).

# FXi-09 Built-In Battery Unit (Ni-MH Rechargeable Battery Pack)

- Charging time: Approximately 10 hours.
- Continuous operation hours: Approximately 8 hours

#### Note

The charging time depends on the operating environment.
 During charging, the balance is not available for use.

### FXi-10 Small Breeze Break


Breeze break shorter than the standard large breeze break.
 (Total height of the small breeze break: 105 mm)


### $FX_{i-1}$ Main Unit Cover

Breeze break provided as standard.

### AX-FXi-31 Main Unit Cover

Main unit protective cover provided as a standard accessory.






### AX-CARAT PAN-WJA Carat Pan (Silver Colored)

 Two silver colored carat pans (the same as the carat pan provided as a standard accessory)


## AX-CARAT PAN-WJA Carat Pan (Black Colored)

 Two black carat pans (the same as the carat pan provided as a standard accessory)



### AX-USB-9P USB Converter

- Adds a COM port to a PC.
- Enables bi-directional communication between the PC and the balance when a USB driver is installed.
- Can use serial communication software such as WinCT on a PC without COM ports.



### AD-1671 Anti-Vibration Table


- A 27-kg table with vibration-absorbing rubber feet, to reduce vibration coming from the floor to enable stable weighing.
- Use the AD-8922A remote controller to avoid weighing errors that can be caused by slight table tilting when operating the scale manually.

## AD-1672 Tabletop Breeze Break

- Protects the balance from wind from sources such as air conditioning or people passing by, thereby reducing balance weighing error.
- The transparent panel assembly consists of antistatic plastic material that protects the balance from static electricity.

### AD-1683 DC Static Eliminator

Used to minimize weighing errors due to static electricity on the material. The AD-1683 is direct-current static eliminator. The ions generated produce no breeze and are effective over a long distance. Therefore, the balance can accurately weight powders, etc. by using the AD-1683.



### AD-1684 Electrostatic Field Meter

This option measures the amount of the static charge on the sample, tare or peripheral equipment and displays the result.
 If those are found to be charged, discharge them using the AD-1683 DC static eliminator.

# AD-1687 Weighing Environment Logger

- A data logger equipped with 4 sensors for temperature, humidity, barometric pressure and vibration that can measure and store environmental data.
- When connected to the RS-232C interface of the balance, the AD-1687 can store environmental data along with weighing data. Therefore, it is possible to store data in an environment where a computer can not be used.
- The stored data can be read to a personal computer using USB.
   As the AD-1687 is recognized as USB memory, special software is not required to read the data.



### AD-1688 Data Logger

- When connected to the RS-232C interface of the balance, the AD-1688 can store the data in an environment where a personal computer can not be used.
- The stored data can be read to a personal computer using USB. As the AD-1688 is recognized as USB memory, special software is not required to read the data.

### AD-1689 Tweezers For Calibration Weight

 This option is used when calibrating the balance using an external weight.

### AD-8121B Printer


- Compact dot-matrix printer
- Statistical function, clock and calendar function, interval print function, graphic print function, dump print mode
- □ 5 x 7 dots, 16 characters per line
- Print paper (AX-PP143, 45 (W) x 50 (L) mm, ø65 mm)
- AC adapter or alkaline battery.

### AD-8920A Remote Display

Connected to the balance using the RS-232C serial interface to display the weight data away from the balance.

### AD-8922A Remote Controller

- Connected to the balance using the RS-232C serial interface to display the weight data and to remotely control the balance.
- Analog output and comparator output available as an option to be installed.





# 20. Specifications

	FZ-1200CT	FZ-700CT	FZ-500CT		
Woighing capacity	1260 ct	760 ct	510 ct		
Weighing capacity	252 g	152 g	102 g		
Maximum display	1260.042 ct	760.042 ct	510.042 ct		
Maximum display	252.008 g	152.008 g	102.008 g		
Minimum weighing value		0.001 ct			
(1 digit)		0.001 g			
Repeatability		0.001 ct			
Standard deviation		0.0005 g			
Linearity		±0.002 ct			
Linearity		±0.001 g			
Stabilization time (typical at FAST)		Approx. 2 second*1			
Sensitivity drift (10 °C - 30 °C / 50 °F - 86 °F )	±2 ppm/°C				
Internal mass		Available			
Clock and calendar function		Available			
Operating environment	5 °C to 40 °C (4	1 °F to 104 °F) 85 % condensation)	%RH or less (No		
Display refresh rate	App	$\frac{1}{10}$ rox. $5^{*2}$ , 10 times/sec	cond		
Counting Minimum unit mass		0.001 g			
mode *3 Number of samples	10	0, 25, 50 or 100 piece	es		
Percent Minimum 100% reference m-ass		0.100 g			
mode *3 Minimum 100% display	0.01%, 0.1%, 1% (I	Depends on the refer	rence mass stored.)		
Interface		RS-232C			
External calibration weight	250 g 200 g 100 g 50 g 25 g	150 g 100 g 50 g 20 g	100 g 50 g 20 g		
Weighing pan diameter	<u> </u>	90 mm			
External dimensions	198(	W) x 294(D) x 315(H)	) mm		
AC adapter	Confirm that the adapter type is correct for the local voltage and power receptacle type				
Power consumption		VA (supplied to the A	•		
Mass of product	Approx. 3.8 kg				

<sup>\*1:</sup> With MID., factory setting of the rate of the response adjustment, the stabilization time is approx. 3 seconds.

<sup>\*2:</sup> Factory setting

<sup>\*3:</sup> To use the counting mode or the percent mode, storing each mode as described in "9–9. Storing Units" is required.

		FZ-1200CT	FZ-700CT	FZ-500CT	
Weighing capacity		1260 ct	760 ct	510 ct	
vveigi iii ig	Сараску	252 g	152 g	102 g	
Maximum	dienlay	1260.042 ct	760.042 ct	510.042 ct	
IVIAXIITIUIT	l display	252.008 g	152.008 g	102.008 g	
Minimum	weighing value	0.001 ct			
(1 digit)		0.001 g			
Repeatab			0.001 ct		
Standard	deviation		0.0005 g		
Linearity			±0.002 ct		
			±0.001 g		
Stabilization (typical at			Approx. 2 second $^{*1}$		
Sensitivity (10 °C - 30	⁄ drift ) °C / 50 °F - 86 °F )		±2 ppm/°C		
Internal m	nass	Unavailable			
Clock and	d calendar function		Unavailable		
Operating	environment	5 °C to 40 °C (41 °F to 104 °F) 85 %RH or less (No condensation)			
Display re	efresh rate	Δnn	rox. $5^{*2}$ , 10 times/sec	rond	
	Minimum unit mass	Αρρ	0.001 g	Oliu	
	Number of samples	1(	0, 25, 50 or 100 piece	<i>5</i> ¢	
Percent	Minimum 100% reference m-ass	-			
mode *3	Minimum 100% display	0.01%, 0.1%, 1% (Depends on the reference mass store			
Interface	iviii iii ii ia ii i i oo /o alopiay	0.01%, 0.1%, 1% (Depends on the reference mass stored.			
torraco		250 g	150 g	100 g	
		200 g	100 g	50 g	
External of	calibration weight	100 g	50 g	20 g	
		50 g	20 g		
		25 g			
	pan diameter	90 mm			
=xiemai c	dimensions		W) x 294(D) x 315(H)		
AC adapt	er	Confirm that the adapter type is correct for the local vo		•	
	nsumption	Approx. 11 VA (supplied to the AC adapter)			
Mass of p	roduct	Approx. 3.4 kg			

<sup>\*1:</sup> With MID., factory setting of the rate of the response adjustment, the stabilization time is approx. 3 seconds.

<sup>\*2:</sup> Factory setting


<sup>\*3:</sup> To use the counting mode or the percent mode, storing each mode as described in "9–9. Storing Units" is required.


# 21. External Dimensions


FZ-1200CT FZ-700CT FZ-500CT FX-1200CT FX-700CT FX-500CT

Unit: : mm

- \*1: When opening maximum
- \*2: Inside dimension
- \*3: Weighing pan diameter
- \*4: Height up to the weighing pan
- \*5: When pulling out the locking handle
- \*6: Inside dimension (from up to the weighing pan)
- \*7: Underhook position (hole center)







# 22. Terms/Index

# 22–1. **Terms**

**Stable value** The weight data when the stabilization indicator appears.

**Environment** Ambient conditions such as vibration, drafts, temperature, static electricity and

magnetic fields which affect the weighing operation.

**Calibration** Adjustment of the balance so that it can weigh accurately.

**Output** To output the weight data using the RS-232C serial interface.

**Zero point** A weighing reference point or the zero display. Usually refers to the value displayed

when nothing is on the weighing pan.

**Digit** Unit of digital resolution. Used for the balance, a unit of minimum weighing value.

**Tare** To cancel the weight of a container which is not included in the weight data.

**Mode** Balance operational function.

**Re-zero** To set the display to zero.

**GLP** Good Laboratory Practice.

**GMP** Good Manufacturing Practice.

**Repeatability** Variation in measured values obtained when the same weight is placed and

removed repetitively. Usually expressed as a standard deviation.

e.g. Standard deviation=1 digit: This means that measured values fall within  $\pm 1$  digit

in the frequency of about 68%.

**Stabilization time** Time required after a sample being placed, until the stabilization indicator illuminates

and the weight data is displayed.

**Sensitivity drift** An affect that a change in temperature causes to the weight data. Expressed as

temperature coefficient.

e.g. Temperature coefficient = 2 ppm/°C : If a load is 200 g and the temperature

changes by 10°C, the value displayed changes by the following value.

0.0002% °C x 10 °C x 200 g = 4 mg

22-2.	Index			
<ul><li>Left, Upper,</li></ul>		22	Breeze break5, 6	. 8
<ul><li>Left, Upper,</li></ul>			Breeze break ring5	•
✓ Left, Upper,	<b>5 5</b>		ЬЕРг Data bit, parity bit	
Left, Upper,			Bubble spirit level	
✓ Left, Lower,			ERL in One-touch calibration	
→ Right, Lower			Calibration	
Right, Lowe			Calibration using an external weight	
RESPONSE FAS	· • • • • • • • • • • • • • • • • • • •		Correcting the internal mass value18 -:	
	Space mark29 - 3		One-touch calibration15,	
IIÓ ON:OFF	ON:OFF key		Calibration report15,	
(ZAL)	CAL key10		Calibration test	
MODE	MODE key		Calibration test using an external weight	17
Q PRINT	PRINT key10		Calibration test report	
◆0← RE-ZERO	RE-ZERO key10		באביי Calibration using an external weight	
1/10d SAMPLE	SAMPLE key10		EANEEL Cancel	
	,	•	[RN Go Cancel Go	43
			[ RN no Cancel No	43
AC adapter	5, 6, 9, 6	9, 70	Carat pan	5
ACAI .		37	[[ ουξ Calibration test	17
AD-1671	Anti-vibration table	66	[[ In Calibration test	15
AD-1672	Tabletop breeze break	66	CE Mark	4
AD-1683	DC static eliminator		Characters for display	47
AD-1684	Electrostatic field meter	66	EL Rad Clock23, 31,	45
AD-1687	Weighing environment logger	3, 67	Cleaning9,	62
AD-1688	Weighing data logger		ELERR Clear	43
AD-1689	Tweezers for calibration weight	67	Clock function45, 69,	70
AD-8121B	Printer5	7, 67	[Lr Cancel	21
AD-8920A	Remote display	67	[Lr [a] Initialization, Perform21,	43
AD-8922A	Remote controller	67	ELr no Initialization, Cancel21,	43
Air buoyancy		8	Command	59
AK code		60	Comparator10,	44
RP Fnc	Application26	6, 40	[ and Condition23, 26, 2	27
ЯР-Ь	Auto print difference24	1, 28	Container11,	37
RPF	Application function26, 39	9, 40	Counting mode33, 34, 37, 69,	70
RP-P	Auto print polarity24	1, 28	Cover	5
Ar-d	Zero after output	24	[P Comparator mode23,	44
AF-E	Auto feed	24	[P Fnc Comparator23,	44
Average		40	EPH Comparator upper limit23,	44
AX-CARAT PA	N-BEX Carat pan, black color	66	[PLa Comparator lower limit23,	44
AX-CARAT PA	N-WEX Carat pan, silver color	66	ErLF Terminator	25
AX-FX <i>i</i> <b>-</b> 31	Cover	66	[5] Internal mass correction18, 19,	26
AX-USB-9P	USB Converter		CSV format	31
685Fnc	Environment & Display23	3, 28	Data bits	
Baud rate		56	Data format29 - 31,	56
6EEP	beep		Data output	47
6EP-	OK buzzer		Data output mode28 - 30,	31
ЬEР⁻	HI buzzer	23	Date	45
PEL <sup>-</sup>	LO buzzer	23	Dimensions69 -	71
bP5	Baud rate	25	Display5,	10

Display refresh	rate	69, 70	Mass of product		69, 70
door		5	Maximum		40
dout	Data output	24, 28, 57	Maximum displa	ay	69, 70
d5	Selection of a sample	: 25	MID.		12
d5 Fnc	Density function	25	Minimum		40
Dust plate		5	Minimum displa	y	69, 70
EMC		4	ñLE	Programmable-unit	25
Environment		63, 72	OL	Overload header	29
Er[d	AK, Error code	25, 60	Operating enviro	onment	69, 70
Error code		64	Output		
FR iL	Fail	55	Pan support		5, 6
FAST		12	Parity		56
FCC		4	PRSS	Password	.53 - 55, 65
Function Table		20, 22, 23	Percent mode	33, 34	, 39, 69, 70
FX <i>i</i> <b>-</b> 02	USB interface	3, 58, 66	PnE	Decimal point	23, 27
FX <i>i</i> -08	Ethernet interface	3, 66	Poff	Automatic display-off	23, 27
FX <i>i</i> -09	Built-in battery	3, 66	P-on	Resume function	
FX <i>i</i> -10	Small breeze break		Power consump	otion	69, 70
FX <i>i</i> -11	Main unit cover	66	Power source '		
GLP			PrE	Data output mode	24, 28
GMP			P5	Function selection mode	
Grounding			PUSE	Data output pause	
•	inal		QT	Percent mode header	
HoLd	Hold function		REU9A		
Hold indicator		•		AD-8920A	
Humidity			Repeatability		
ıd	ID number setting		RESPONSE		
ID label			rnli	Dispaly at start	•
ID number			RS-232C		
inFo	GLP output		RsCom	Command mode for WinCT	
Initializing		•	RsKey	Data output mode for WinCT	58
ınt	Interval time		RsWeight	Data output mode for WinCT	
Interface			Sensitivity drift		
Interval mode	in		5- ıd	ID number output	
Key			5 iF	Serial interface	
•		10. 22	SLOW		
MODE Key		,	5Pd	Display refresh rate	
,	(ey		ST	Stable header	
			Stability		
	ey		•	j	
	ey		Stable value		
Key Operation			SERF	Statistic display	
Keys				on	
Ld in	Liquid density input		Static electricity		
Leveling foot	Elquid deriolty input		•		
Linearity			5t-b	Stability detection	
Lo	Sample mass error	,	5-Ed	Time & Date output	
Loc na	Password input		Stop bit	Time & Date output	
LocFnc	Password function		Stream mode		
	ance	•	Sum		
way iono subsu	ы ю		Juin		

SurE	Sure53 - 55	<u>LL</u>	Tael (Preset in the factory)33, 34
Symbols	10	EoL	Tola (India)33, 34
Tare	11	MES	Messghal33, 34
Temperature	8	Lb	Pound33, 34
Terminator	29 - 31	L OZ	Pound Ounce33, 34
Time	45	MLT	Programmable unit33, 34, 52
Erc	Zero tracking23, 27	Un ıE	Unit mode25
E-UP	Timeout25	URL	63
E YPE	Data format25	US	Unstable header29, 30
Underhook	8, 51	USB	3
Unit		Weighing capa	city69, 70
g	Gram33, 34	Weighing mode	e11
PE <b>S</b>	Counting mode33, 34, 37	Weighing opera	ation11
%	Percent mode33, 34, 39		5, 6, 11, 37
$\Box Z$	Ounce (Avoir)33, 34	Weight or Mass	s13, 16 - 19
□Zt	Troy Ounce33, 34	WinCT	Communication Tools58
Ľť	Metric carat33, 34	WT	Counting mode header29
m <u>o</u> m	Momme33, 34	WT	Stable header, DP format30
dint.	Pennyweight33, 34	Zero point	72
	1 011119 4401911111111111111111111111111	Zoro pon k	

# **MEMO**



### A&D Company, Limited

3-23-14 Higashi-Ikebukuro, Toshima-ku, Tokyo 170-0013, JAPAN Telephone: [81] (3) 5391-6132 Fax: [81] (3) 5391-6148

#### **A&D ENGINEERING, INC.**

1756 Automation Parkway, San Jose, California 95131, U.S.A. Telephone: [1] (408) 263-5333 Fax: [1] (408)263-0119

### **A&D INSTRUMENTS LIMITED**

Unit 24/26 Blacklands Way, Abingdon Business Park, Abingdon, Oxfordshire OX14 1DY United Kingdom Telephone: [44] (1235) 550420 Fax: [44] (1235) 550485

#### A&D AUSTRALASIA PTY LTD

32 Dew Street, Thebarton, South Australia 5031, AUSTRALIA Telephone: [61] (8) 8301-8100 Fax: [61] (8) 8352-7409

#### **A&D KOREA Limited**

한국에이.엔.디(주)

대한민국 서울시 영등포구 여의도동 36-2 맨하탄 빌딩 8층 우편 번호 150-749 (Manhattan Building 8th Floor, 36-2 Yoido-dong, Youngdeungpo-gu, Seoul, 150-749 Korea)

전화: [82] (2) 780-4101 팩스: [82] (2) 782-4280

#### OOO A&D RUS

ООО "ЭЙ энд ДИ РУС"

121357, Российская Федерация, г.Москва, ул. Верейская, дом 17 (Business-Center "Vereyskaya Plaza-2" 121357, Russian Federation, Moscow, Vereyskaya Street 17)

тел.: [7] (495) 937-33-44 факс: [7] (495) 937-55-66

### **A&D INSTRUMENTS INDIA PRIVATE LIMITED**

ऐ&डी इन्स्ट्रयमेन्ट्स इण्डिया प्रा0 लिमिटेड

509, उद्योग विहार , फेस -5, गुड़गांव - 122016, हरियाणा , भारत

( 509, Udyog Vihar, Phase-V, Gurgaon - 122 016, Haryana, India ) फोन : 91-124-4715555 फैक्स : 91-124-4715599