GX-M/GF-M

Quick Start Guide
GX-M/GF-M Series
High capacity precision balances

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© 2020 A&D Company, Limited. All rights reserved.

No part of this publication may be reproduced, transmitted, transcribed, or translated into any language in any form by any means without the written permission of A&D Company, Limited.

The contents of this publication and the specifications of the instrument covered by this publication are subject to change for improvement without notice.

Quick Start Guide – GX-M/GF-M Series	EN
High capacity precision balances	4
Guide de démarrage rapide – Séries GX-M/GF-M	FR
Balance de précision grande portée	14
Schnellstartanleitung – Serie GX-M/GF-M	DE
Hochlastige Präzisionswaage.....	24
Guía de inicio rápido – Series GX-M/GF-M	ES
Balanzas de precisión de alta capacidad	34
Guida introduttiva – Serie GX-M/GF-M	IT
Bilance di grande portata e precision	44
快速入門指南- GX-M/GF-M 系列	CT
大秤量精密天平.....	54

1.Introduction

This manual describes the basic functions and operations of the balance.

For more details, refer to each manual listed in "1-1. List of Functions Explained in Separate Manuals".

1-1 List of Functions Explained in Separate Manuals

The functions and operating procedures of the balance are explained in detail in other documents besides this manual.

The following three manuals can be downloaded from our website <https://www.aandd.jp>.

1. [Communication Manual]

Using the communication function of the balance, a personal computer (PC) or printer

This is a supplementary manual for connecting to peripheral devices.

2. [GX-M series/GF-M series Instruction Manual]

This is an instruction manual for understanding and fully utilizing the precision electronic balance GX-M/GF-M series.


3. [Flow Rate Display (FRD) Function]

This is a supplementary manual for using the flow rate measurement function

(FRD: Flow rate display).

2. Installation And Precautions

The balance is a precision instrument. Unpack it carefully. Keep the packing material to be used for transporting the balance in the future. The packing contents depend on the balance model. See the illustrations to confirm that everything is contained.

Breeze break provided for
GX-8202M, GX-8202MD, GX-10202M
GF-8202M, GF-8202MD, GF-10202M


Note

- Please confirm that the AC adapter type is correct for your local voltage and receptacle type.
- Please use the specified dedicated AC adapter for the balance.
- Do not use the included AC adapter for models that are not considered compatible with the AC adapters.
- If you use the wrong AC adapter, the balance and other equipment may not operate properly.

2-1 Precautions

To get the optimum performance from the balance and acquire accurate weighing data, note the following:

- Install the balance in an environment where the temperature and humidity are not excessive. The best operating temperature is about $20^{\circ}\text{C}\pm 2^{\circ}\text{C}$ at about 45 ~ 60%RH relative humidity.
- Install the balance where it is free of dust.
- The weighing table should be solid, free from vibration and drafts, and as level as possible. (An anti-vibration table or stone table is ideal)
- Install the balance in a stable place avoiding vibration and shock. Corners of rooms on the first floor are best, as they are less prone to vibration.
- Install the balance where it is not affected by heaters or air conditioners.
- Install the balance where it is not exposed to direct sunlight.
- Install the balance away from equipment which produces magnetic fields.
- Level the balance by adjusting the leveling feet and confirm it using the bubble spirit level.
- Warm up the balance for at least 30 minutes. Plug in the AC adapter as usual.
- Calibrate the balance before use or after moving it to another location. Refer to the separate instruction manual "9. Calibration".


Caution

Do not install the balance where flammable or corrosive gas is present.


3.Display Symbols And Key Operation

Display symbols

- Number of statistical data (Statistical calculation mode)
- Displays the weight data relative to the weighing capacity, in percentage, in the weighing mode (**Capacity indicator**)

Blinking display contents


Key operation

Key operations affect how the balance functions. Normal key operation during measurement is “Press and release the key immediately” or “Press the key”. Please do not “Press and hold the key (for 2 seconds)” unless required.

Press the key
(Press and release
the key immediately)

Press and hold
the key for 2
seconds

Key	When pressed and released	When pressed and held (for 2 seconds)
	Turns the display ON:OFF . The standby indicator is displayed when the display is turned off. The weighing mode is enabled when the display is turned on. When password function is enabled, password input display will be displayed. Refer to the separate instruction manual "21-2 How to Input The Password At The Start Of Weighing" This ON:OFF key is available anytime. Pressing the ON:OFF key during operation will interrupt operation and turn the display OFF. *	
	In the weighing mode, turns the digit for readability on and off. In the counting or percent mode, enters the sample storing mode.	<ul style="list-style-type: none"> - Enters the function table mode. Please refer to the separate instruction manual "11. Function Table". - Runs the repeatability check function when pressed and held for another 2 seconds after the function table menu is displayed. Please refer to the separate instruction manual "22. Repeatability Check Function". (GX-M series only)
	Switches the weighing units stored in the function table.	Displays other items of the calibration menu. Please refer to the separate instruction manual "8-2 Self Check Function/Automatic Setting Of Minimum Weighing Value by ECL".
	Performs calibration of the balance using the internal mass. (GX-M series)	Displays other items of the calibration menu.
	Stores the weighing data in memory or outputs to a printer or personal computer depending on the function table settings. (Factory setting = output)	Enters mode to change the unit mass registration number in counting mode. By changing the function table: <ul style="list-style-type: none"> - Outputs "Title block" and "End block" for GLP, GMP report. - Displays the data memory menu. - Enters mode for reading density number in flow measurement.
	Sets the display to zero.	

* When the "Gross net tare function" is selected, the display is turned off by pressing and holding (for 2 seconds). Please refer to the separate instruction manual "16. Gross Net Tare Function".

4. Calibration

Since the balance's resolution is high, weighing values may change due to gravity and daily environmental changes.

4-1 Automatic Calibration (GX-M Series Only)

The GX-M series automatically adjusts the balance sensitivity using the built-in weight according to the temperature change of the operating environment or the set time and interval time. For details, refer to "9-1 Auto calibration" in the separate instruction manual.

4-2 One-Touch Calibration (GX-M Series Only)

This function calibrates the balance using the internal mass.




1. Connect the AC adapter and warm up the balance for at least 30 minutes with nothing on the weighing pan.
2. Press the **CAL** key. The balance displays **[CAL in]**.
3. The balance performs calibration using the internal mass. Do not allow vibration or drafts to affect the balance.
4. After the calibration, If the "GLP output (inFd)" parameter of the function table is set, the balance outputs a sensitivity adjustment execution record.
5. The balance will automatically return to weighing mode after calibration.









4-3 Calibration Using An External Weight

This function calibrates the balance using an external weight.

1. Connect the AC adapter and warm up the balance for at least 30 minutes with nothing on the weighing pan.
2. Press and hold the **CAL** key for 2 seconds until **[CAL OUT]** is displayed, then release the key.
3. Make sure that nothing is on the weighing pan and press the **PRINT** key to weigh the zero point. Do not apply vibration, etc.
4. Place the external weight on the weighing pan and press the **PRINT** key. Do not apply vibration etc.
5. Remove the external weight from the weighing pan.
6. After calibration, if GLP output is set, "sensitivity adjustment execution record" is output or stored in data memory.
7. The display automatically returns to weighing display.
8. Place the external weight again and check that the set value is ± 2 digits.
If it is out of range, pay attention to the surrounding environment and start from "1".

4-4 Error Codes

Display	Error code	Description
	EC, E11	Stability error The balance cannot stabilize due to an environmental problem. Check around the pan. Prevent vibration, drafts, temperature changes, static electricity and magnetic fields, from influencing the balance. Refer to the separate instruction manual "3-2 During Use". To return to the weighing mode, press the  key.
		Out of the setting range The data to be stored is out of the setting range.
		Malfunction of the internal memory element of the balance If this error continues to be displayed, repair is necessary. Please contact your dealer.
	EC, E16	Internal mass error Applying the internal mass does not yield a change in the mass value as specified. Confirm that there is nothing on the pan and perform the weighing operation from the beginning again.
	EC, E17	Internal mass error The internal mass application mechanism does not function properly. Perform the weighing operation from the beginning again.
		Abnormality in the internal memory data of the balance If this error continues to be displayed, repair is necessary. Please contact your dealer.
		Abnormality in the internal memory data of the balance If this error continues to be displayed, repair is necessary. Please contact your dealer.
	EC, E20	Calibration weight error (Positive value) The calibration weight is too heavy. Confirm the calibration mass value. Press the  key to return to the weighing mode.
	EC, E21	Calibration weight error (Negative value) The calibration weight is too light. Confirm the calibration mass value. Press the  key to return to the weighing mode.
		Overload error A sample beyond the balance weighing capacity has been placed on the pan. Remove the sample from the pan.
		Weighing pan error The mass value is too light. Confirm that the weighing pan is properly installed and calibrate the balance.

Display	Error code	Description
		<p>Sample mass error</p> <p>The balance cannot store the sample for the counting mode or for the percent mode because it is too light. Use a larger sample.</p>
  		<p>Unit mass error</p> <p>The sample unit mass for the counting mode is too light. Storing and using it for counting will cause a counting error. Add samples to reach the specified number and press the PRINT key. Pressing the PRINT key without adding samples will shift the balance to the counting mode. But, for accurate counting, be sure to add samples.</p>
		<p>Clock battery error</p> <p>The clock backup battery has been depleted. Press any key and set the time and date. The clock and calendar function works normally as long as the AC adapter is connected to the balance. If this error appears frequently, contact the local A&D dealer.</p>
		<p>Power supply voltage fault</p> <p>The voltage supplied from the AC adapter is abnormal. Please check if the problem is the AC adapter supplied with the balance.</p>
 		<p>ECL repeatability</p> <p>With the self-check function, the standard deviation (SD) of repeatability due to electronically controlled load (ECL) exceeded 50 digits. Please revise the installation environment of the balance.</p> <ul style="list-style-type: none"> <input type="checkbox"/>  This is displayed when repeatability is displayed by ECL. <input type="checkbox"/>  This is displayed when the minimum weighing value (reference value) by ECL is displayed. <p>Refer to the separate instruction manual "8-2 Self Check Function / Automatic Setting Of Minimum Weighing Value by ELC.</p>
		<p>Full memory</p> <p>The maximum number or stored weighing values has been reached. In order to store more weighing values, it is necessary to delete the data. Refer to the separate instruction manual "13. Data Memory".</p>
		<p>Full memory</p> <p>The stored calibration and calibration test history has reached 50 instances. If more is stored, the old history will be deleted. Refer to the separate instruction manual "13. Data Memory"</p>

Display	Error code	Description
	EC, E00	Communications error A protocol error occurred in communications. Confirm the format, baud rate and parity.
	EC, E01	Undefined command error An undefined command was received. Confirm the command.
	EC, E02	Not ready A received command cannot be processed. Example: The balance received a "Q" command, but not in the weighing mode. The balance received a "Q" command while processing a RE-ZERO command. Adjust the delay time to transmit a command.
	EC, E03	Timeout error If the timeout parameter is set to $t-UP$ 1^{SEC} , the balance did not receive the next character of a command within the time limit of one second. Confirm the communication.
	EC, E04	Excess characters error The balance received excessive characters in a command. Confirm the command.
	EC, E06	Format error A command includes incorrect data. Example: The data is numerically incorrect. Confirm the command.
	EC, E07	Parameter setting error The received data exceeds the range that the balance can accept. Confirm the parameter range of the command.

4-5 Other Display

Advice

When this indicator (◀) blinks, automatic self calibration is waiting. If the balance is not used for several minutes with this indicator blinking, the balance automatically performs calibration using the internal mass. The blinking period depends on the operating environment.

The balance can be used when this indicator is blinking, but we recommend that you perform calibration before weighing.

5. Specifications

	GX-8202M GF-8202M	GX-8202MD GF-8202MD	GX-10202M GF-10202M
Weighing capacity	8.2kg		10.2kg
Maximum display	8.20084kg	8.2008kg	10.20084kg
		2.20009kg *1	
Readability	0.01g	0.1g	0.01g
		0.01g	
Repeatability (Standard deviation)	0.01g	0.05g	0.01g
		0.01g	
Linearity	±0.03g	±0.1g	±0.03g
		±0.02g	
Accuracy after calibration with internal mass *2	±0.15g	±0.3g	±0.15g

	GX-12001M GF-12001M	GX-22001M GF-22001M	GX-32001M GF-32001M	GX-32001MD GF-32001MD
Weighing capacity	12.2kg	22.2kg	32.2kg	
Maximum display	12.2084kg	22.2084kg	32.2084kg	32.208kg
				6.2009kg*1
Readability	0.1g			1g
				0.1g
Repeatability (Standard deviation)	0.1g			0.5g
				0.1g
Linearity	±0.2g			±1g
				±0.2g
Accuracy after calibration with internal mass *2	±1.0g		±1.5g	±3g

*1 The balance allows weighing using a precision range, even with a heavy tare placed on the pan.
(Smart range function)

*2 The operating environment does not include excessive change of ambient temperature, humidity,
vibration, drafts, magnetic fields and static electricity.

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1756 Automation Parkway, San Jose, California 95131, U.S.A.
A&D AUSTRALASIA PTY LTD
A&D KOREA Limited