

Rigid-body Pendulum type Physical Properties Testing Instrument

ND RPT-3000






RPT-3000W Highly accurate measured

Installation example of small UV irradiation device (optional)


UV irradiation from both sides of the knife edge.




At an angle.



How to read measurement data



How to Read Physical Properties Data



Introduction

The Rigid-body Pendulum Type Physical Properties Testing Instrument is an instrument that dynamically measures the curing process of a substance and the physical properties of the surface. The behavior of the various raw materials used to form a solidified coating film has a great influence on the characteristics of the

subsequent film.


Also, the mechanical properties of the finished coating film are a major concern.

RPT

The solution-like and powder-like substances on the base material solidify undervarious conditions (temperature, humidity, light, time, etc.) to form a coating film. Even if the solidified material looks fine. if those conditions are not examined properly, then the characteristics of the solidified material cannot be understood properly, which may cause problems later.

If the curing conditions are decided, then it may become necessary to select appropriate materials and compositions for those curing conditions or may involve research and development of new materials. In addition, it can be used as quality control before it goes on the market by measuring the surface characteristics of the finished coating film. RPT-3000W is an indispensable item for such people. ISO 12013-1 was established for the examination of curing conditions, and ISO 12013-2 was established for surface physical characteristics, which can be measured with RPT-3000W.






Example of surface physical property measurement of PS thin film (50, 100, 200 nm) (Measurement data for example was provided by Dr. Inuzuka of Tanaka Laboratory Kyushu University Graduate School)

Rigid-body Pendulum type Physical Properties Testing Instrument

ment of drying/curing behavior and physical properties of substances.

Measurement Mode (curing, physical characteristics)

In a curing measurement, a knife edge is placed in a sample on the base material which displaces the knife attached pendulum very slightly. In surface physical characteristics measurement, a cylinder edge is placed on the measurement surface (see the figure below). By releasing the applied displacement, the pendulum vibrates freely while being influenced by the sample. The cycle period and logarithmic decrement is found.






Rigid-body pendulum



Heating/cooling block (CHB-100, CHB-200)




Substrate, edge





Coating Tool (CT)

Measurement Examples using RPT


Measurement Examples of Curing Behavior








Measurement Examples of Physical Characteristics



Effect of body soap on artificial skin

[C]



RPT can be used to evaluate materials in a wide range of fields

Paint Adhesives	 Evaluation of curing temperature and curing time. Evaluation of curing/physical properties of solidified materials using curing agents. Design data for painting/coating production lines. Quality evaluation of coating films and thin films (50 nm). Evaluation of adhesiveness, distortion, etc. Early prediction of weather resistance deterioration. Others
Cosmetics Pharmaceuticals	 Evaluation of dryness and surface properties of nail polish and mascara. Evaluation of dryness, adhesiveness, and elasticity of packs and poultices. Evaluation of lubricity and cleansing properties of hair, etc. Evaluation of dryness of contact lenses. Others
Food	 Evaluation of gelling properties of gelatin, agar, etc. Quantitative evaluation of sensory performance of foods, etc. Others
Fibers	 Evaluation of physical properties of fibers. Evaluation of the texture of the cloth. Others

Plastics	 Evaluation of surface physical properties. Evaluation of film physical properties. Hard coat curing, surface and internal physical characteristics evaluation. Others
Printing	 Evaluation of ink physical properties on printed matter. Evaluation of ink drying property. Evaluation of transferability to rollers. Others
Electrical Electronics	 Evaluation of physical properties of battery materials (electrode film, spacers, etc.). Evaluation of curability and physical properties of optical fibers and optical filters. Evaluation of curing characteristics of conductive paste. Evaluation of melting and solidification characteristics of solder. Curing of printed circuit board and evaluation of physical properties. Others
Others	• Evaluation of physical properties of concrete, asphalt, etc.

Main Specifications

Amplitude displacement detection: Non-contact eddy current displacement sensor Maximum amplitude angle: ±0.57 degrees

Angle resolution: 1.75 e⁻⁵ degree

Vibration cycle: 0.050 to 2.000 seconds

Logarithmic decrement: 0.001 to 3.0 Measurement temperature range: -80 to +400°C

Cooling method: Liquid nitrogen

Optional UV irradiation device, etc.

External dimensions/weight: Main unit 300 x 220 x 525 mm / 15 kg

(WxDxH) Control unit 410 x 350 x 135 mm / 12 kg Safety device: Overheat prevention device

Warning lamp (ON at 50°C or higher)

Power supply AC 100 V, 550 VA

Software: OS Windows 10 Pro

Basic application MSAT0001V2 Application MSAT0010V2 The Rigid Pendulum Type Physical Properties Testing Instrument (corresponding to RPT-3000W) owned by A&D Co., Ltd., has been adopted as the international standard in ISO 12013-1 and ISO 12013-2. (October 2012)

ISO 12013-1: How to measure paint cross linkage and mesh formation temperature. ISO 12013-2: Method for measuring thermal properties (Tg, physical properties, etc.) of paint.



• Please read the instruction manual carefully and use the equipment correctly.



Discover Precision

A&D Company, Ltd.

3-23-14 Higashi-Ikebukuro, Toshima-Ku, Tokyo, 170-0013, Japan Tel: +81 3-5391-6132 Fax: +81 3-5391-1566 http://www.aandd.jp

A&D Engineering, Inc. 1756 Automation Parkway, San Jose, CA 95131, U.S.A. Tel: +1 408-263-5333 Fax: +1 408-263-0119

A&D Korea Ltd.

8F Manhattan Bldg., 33, Gukjegeumyung-ro 6-gil, Yeongdeungpo-gu, Seoul, 07331, Korea Tel: +82 2-780-4101 Fax: +82 2-782-4280

A&D Technology Trading Co.,Ltd.

32CD,World Plaza,No.855 South Pudong Road, China(Shanghai) Pilot Free Trade Zone,200120,China Tel: +86 21-3393-2340 Fax: +86 21-3393-2347

A&D Scientech Taiwan Ltd.

4F No.5 Ching Tao East Road, Taipei, Taiwan, R.O.C. Tel: +886 2-2322-4722 Fax: +886 2-2392-1794

NOTE: Appearance and specifications are subject to change without notice for improvement. • The contents of the catalog are current as of December 2019. * RPT3000W-ADJC-14-ZW9-19c01GP