

Training Material for Balances (2)

Prepared by R&D Division 1

A&D Company, Limited

Contents

- Measuring Repeatability, Linearity & Sensitivity Drift - P.3
 - Specifications of GX-2000 ----- P.4
 - Repeatability ----- P.5
 - Linearity ----- P.10
 - Sensitivity Drift ----- P.12

Measuring *Repeatability*, *Linearity* & *Sensitivity Drift*

Specifications of GX-2000

Weighing Capacity	2100g
Minimum Weighing Value (1 digit)	0.01g
Repeatability (Standard Deviation)	0.01g
Linearity	±0.02g
Sensitivity Drift (10 ~ 30°C)	±2ppm/°C (Automatic self calibration is not used)





Repeatability

- **Repeatability** indicates the variation in measurement when the same mass is loaded several times. Repeatability is expressed using standard deviation (σ : sigma).
- As the repeatability of GX-2000 is 0.01g, when a 2kg mass is loaded repeatedly, the results will be*
 - 1999.99g ~ 2000.01g with the probability of 68% (Mean $\pm 1\sigma$)
 - 1999.98g ~ 2000.02g with the probability of 95% (Mean $\pm 2\sigma$)
 - 1999.97g ~ 2000.03g with the probability of 99.7% (Mean $\pm 3\sigma$)

* supposing that the mean value is 2000.00g

Repeatability

Normal Distribution



Repeatability

Calculating Repeatability (σ)

Formula: Standard Deviation $\sigma_{n-1} = \sqrt{\frac{\sum (\text{Measured Value} - \text{Mean})^2}{\text{Number of Measurements} - 1}}$

Ex) We measured a 2kg mass five times with GX-2000 (minimum weighing value 0.01g). The results were as follows:

1st	2000.01g
2nd	2000.00g
3rd	1999.98g
4th	1999.99g
5th	2000.02g

Repeatability

	Measured Value	Minimum Weighing Value (digit)
1st	2000.01g	1 dig
2nd	2000.00g	0 dig
3rd	1999.98g	-2 dig
4th	1999.99g	-1 dig
5th	2000.02g	2 dig
Mean	2000.00g	0 dig

As the minimum weighing value of GX-2000 is 0.01g,

1 digit = 0.01g

$$\sigma_{n-1} = \sqrt{\frac{\sum (\text{Measured Value} - \text{Mean})^2}{\text{Number of Measurements} - 1}}$$

Measured Value - Mean	(Measured Value - Mean) ²
0.01g	0.0001
0.00g	0.0000
-0.02g	0.0004
-0.01g	0.0001
0.02g	0.0004

$$\sum (\text{Measured Value} - \text{Mean})^2 = 0.0010$$

$$\frac{\sum (\text{Measured Value} - \text{Mean})^2}{\text{Number of Measurements} - 1} = 0.00025$$

$$\sigma_{n-1} = \sqrt{\frac{\sum (\text{Measured Value} - \text{Mean})^2}{\text{Number of Measurements} - 1}} = 0.01588 = 0.016g$$

Repeatability

Alternatively,

Measured Value - Mean	(Measured Value - Mean) ²
1 dig	1
0 dig	0
-2 dig	4
-1 dig	1
2 dig	2

$$\sum (\text{Measured Value} - \text{Mean})^2 = 10$$

$$\frac{\sum (\text{Measured Value} - \text{Mean})^2}{\text{Number of Measurements} - 1} = 2.5$$


$$\sigma_{n-1} = \sqrt{\frac{\sum (\text{Measured Value} - \text{Mean})^2}{\text{Number of Measurements} - 1}} = 1.588 \doteq 1.6 \text{ dig} = 0.016 \text{ g}$$

Since the repeatability of GX-2000 is 0.01g, the performance of this balance would be **unsatisfactory**.


Linearity

- **Linearity** refers to the deviation in the measured value from the ideal (linear) curve between 0g and the maximum capacity.
- It is necessary to correct the slope of the measurement curve by conducting a calibration first, because the weighing results will vary depending on the acceleration of gravity at the point of use or can be influenced by temperature changes.

Linearity



⇒
CAL



- When calibrated with a 2000g weight, the measurement curve will be corrected by -0.04g at the point of 2000g.
- Likewise, it will be corrected by -0.01g , -0.02g , and -0.03g at 500g, 1000g, and 1500g respectively. The maximum linearity error after the calibration will thus be -0.02g (at 1000g).
- Since the specification of linearity for GX-2000 is $\pm 0.02\text{g}$, the performance is judged to be satisfactory in this case.

Sensitivity Drift

- **Sensitivity drift** is the increase or decrease in the measured value that occurs when the same mass is weighed at different ambient temperatures.
- The rate of variation is expressed in “ppm (Parts Per Million).”
- For example, the sensitivity drift of GX-2000 is $\pm 2\text{ppm}/^{\circ}\text{C}$. This means that the measured value can vary by up to “ $\pm 2/1,000,000$ ” with a 1°C change in temperature.

Sensitivity Drift

Ex 1) We calibrated GX-2000 so that it would display 2000.00g at 20°C. If the ambient temperature later became 30°C, what would be the maximum error of measurement caused by sensitivity drift?

Measured Value (Referential Value) × (ppm / °C) × Variation in Temperature

$$\begin{aligned} & 2000.00(\text{g}) \times \{\pm 2(\text{ppm} / ^\circ\text{C})\} \times \{30(^{\circ}\text{C}) - 20(^{\circ}\text{C})\} \\ &= 2000.00 \times \frac{\pm 2}{1000000} \times 10 \\ &= \pm 0.04(\text{g}) \end{aligned}$$

This indicates that there can be a variation of up to $\pm 0.04\text{g}$ (1999.96 ~ 2000.04g) from the initial value of 2000.00g when the temperature changes by 10°C.

Sensitivity Drift

Ex 2) The following results were obtained using GX-2000. What is the sensitivity drift (ppm/°C)?

Weight	Ambient Temperature	Measured Value
2000.00g	20°C	2000.00g
	30°C	1999.97g

$$\left(\frac{\text{Result B (after change)} - \text{Result A (before change)}}{\text{Result A (before change)}} \right) \times 1000000 = (\text{ppm})$$

$$\frac{\text{ppm}}{\text{Variation in Temperature}} = (\text{ppm}/^{\circ}\text{C})$$

$$\frac{1999.97 - 2000.00}{2000.00} \times 1000000 = -15(\text{ppm})$$

$$\frac{-15}{30 - 20} = -1.5(\text{ppm}/^{\circ}\text{C})$$

This satisfies the specification of sensitivity drift for GX-2000, which is $\pm 2\text{ppm}/^{\circ}\text{C}$.