

DIGITAL BLOOD PRESSURE MONITOR

UA-704

Clinically
Validated

IHB



Palm Top Monitor

- Innovative “Palm Top” Style
- Irregular Heart Beat (IHB) Indicator
- SlimFit Comfortable Cuff
- Uses Single AA Battery

CE 0123

A&D
A&D Medical

Technology You Can Trust

<http://www.aandd.jp>

DIGITAL BLOOD PRESSURE MONITOR UA-704

Innovative “Palm Top” Style

An innovative “Palm Top” style makes it easy to monitor your blood pressure. The ergonomically designed monitor fits in your palm and you can touch the START button without any constriction.

Irregular Heart Beat (IHB) Indicator

This lights up if IHB is detected during the measurement. Helps users to find an irregular heartbeat, which may cause cardiovascular disease.

(Note: This is for warning purposes only. A&D recommends contacting a physician if the user sees this indicator frequently)

SlimFit Comfortable Cuff

A new innovative edge free internal bladder “EFB” promises even pressure on the arm during the measurement. This provides the user with a higher level of comfort while the arm is under pressure.


Uses Single AA Battery

The monitor is designed to conserve energy, promising a long unit service life to keep your measurements going.

(Battery life: 2,000 measurements)

WHO blood pressure classification

Standards for assessment of high blood pressure, without regard to age, have been established by the World Health Organization (WHO), as shown in the chart below.



A&D
A&D Medical

Technology You Can Trust

<http://www.aandd.jp>

*UA704CA-ADCC-08-PR5-143**
REF*CARD-UA704