

Notification of Change for ABPM TM-2440/2441

March, 2019

1) Summary

Please be informed that we will change the specifications of the A&D TM-244x Series (Ambulatory blood pressure monitors) as follows:


- The color of ABPM cuff, cuff cover and cuff cloth will be changed to dark navy from the existing blue-gray color.
- The current rechargeable battery (AS-3UTG-2BP) and adaptor (AS-NC-M58) will no longer be used.

2) Reasons for these Changes

Cuff color: Production process and change of raw materials

Battery and adaptor: End of production

3) Appearance of the cuff before and after change



4) Timing of change

March 2019

5) Applicable models

TM-2440, TM-2441

Option Items

Order code	Name	Description
TM-CF502A	Extra large cuff for left arm	Arm circumference 36 to 50 cm 14.2" to 19.7"
TM-CF402A	Large cuff for left arm	Arm circumference 28 to 38 cm 11.0" to 15.0"
TM-CF302A	Adult cuff for left arm	Arm circumference 20 to 31 cm 7.8" to 12.2"

TM-CF202A	Small cuff for left arm	Arm circumference 15 to 22 cm	5.9" to 8.7"
TM-CF802A	Adult cuff for right arm	Arm circumference 20 to 31 cm	7.8" to 12.2"
AX-133024503-S	Extra large cuff cover	for left arm	10 sheets
AX-133024663-S	Large cuff cover	for left arm	10 sheets
AX-133024500-S	Adult cuff cover	for left arm	10 sheets
AX-133024667-S	Small cuff cover	for left arm	10 sheets
AX-133024353-S	Adult cuff cover	for right arm	10 sheets
AX-133025103-S	Extra large cuff cloth	for left arm	2 sheets
AX-133025102-S	Large cuff cloth	for left arm	2 sheets
AX-133024487-S	Adult cuff cloth	for left arm	2 sheets
AX-133025101-S	Small cuff cloth	for left arm	2 sheets
AX-133025104-S	Adult cuff cloth	for right arm	2 sheets

END