

Bedside Scale

Model: AD-6121B

The AD-6121B is a lightweight, portable scale ideal for bedside use.

General Specifications

- Can display measurements in kg
- Lightweight 3 kg design with handle for portability and convenient storage
- Use at the bedside and store underneath
- Adjustable indicator
- Auto power off function
- Supports AA battery operation or AC adapter connection



Level adjustment

Lightweight & compact size



Technical Specifications

Capacity	150 kg
Graduation	0.1 kg
Measurement Range	2 kg to 150 kg
Power Source	4 Batteries LR6 (AA) / AC Adapter
Battery Life	Approx. 1,000 measurements (Battery use)
Interface	RS 232C 1ch (3Pin)
Dimensions (W) x (D) x (H)	315 x 314 x 57 mm
Weight	3 kg (Unit without batteries)
Functions and Properties	TARE, Auto power-off

Optional Accessories

Product Code	Option Name
AX-KO2238	Communication cable
AD-6121-01	Carry case



Communication cable



Carry case

What is In The Box

- 1 Instruction manual
- 1 Warranty

- Please purchase an AC adapter that complies with the standards of your region.
- Ensure that the product complies with the regulatory requirements of the country where it will be used before use.
- Specifications may differ in each country. Please contact an A&D subsidiary for more information.

Baby Scale “Baby One”

Model: AD-6020-12K

The AD-6020-12K is a baby scale with a removable basket for easy cleaning.

General Specifications

- Can display measurements in kg
- Ensure hygiene with easy-to-clean removable basket
- Large 25 mm digital display
- Supports D battery operation or AC adapter connection
- Record and recall last measurement values with a single button
- Tare function



Memory
Memory Function

888
LCD Display


OFF
Auto Power Off

Battery
Normal Battery

Height measurement

The basket has a scale and a cursor for easy height measurement. The cursor for height measurement can be stored in the lower part of the main unit.

Cursor for height measurement



Technical Specifications

Capacity	12 kg
Graduation	10 g < 10 kg > 20 g
Measurement Range	0.1 kg to 12 kg
Height Measurement Range	45 cm to 57 cm
Power Source	4 Batteries LR20 (D) / AC Adapter
Dimensions (W) x (D) x (H)	628 x 343 x 245 mm
Weight	5 kg (Unit without batteries)
Functions and Properties	TARE, Height measurement, Safety lock



What is In The Box

- 1 Height measurement cursor
- 1 Instruction manual
- 1 Warranty

- Please purchase an AC adapter that complies with the standards of your region.
- Ensure that the product complies with the regulatory requirements of the country where it will be used before use.
- Specifications may differ in each country. Please contact an A&D subsidiary for more information.