


IPX5
Water
proof

Pocket pH Meter AP-21

A&D
A&D Weighing

Performs in various fields including
food and beverages, agriculture, water quality, education, and industry



Measurement range: pH 0.0 to 12.0

Product Details



- 1- to 3-point pH calibration
- HOLD function
- Sensor sensitivity display function
- Automatic Temperature Compensation (ATC) function
: Compensation within the range of 0 °C to 50 °C

Specification	AP-21
Measurement method	Glass electrode method
Measurement range	pH: 0.0 to 12.0 pH (Display range: 0.0 to 14.0 pH), Temperature: 0.0 to 50.0 °C
Resolution	pH: 0.1 pH, Temperature: 0.1 °C
Measurement accuracy	pH: ±0.2 pH, Temperature: ±1.0 °C
Operating environment	Temperature: 0 to 50 °C, Humidity: 80 %RH or less (no condensation)
Calibration	1- to 3-point calibration
Automatic temperature compensation	Yes (0 °C to 50 °C)
Auto power off	Approx. 20 minutes
Water resistance	IPX5-equivalent
Power source / Continuous measurement time	4 x LR44 alkaline button cell batteries / Approx. 200 hours
Dimensions / Weight	Approx. Φ38 x 179 (L) mm / Approx. 88 g
Accessories	Battery (for testing), hand strap, instruction manual

A&D
A&D Company, Ltd.

Discover Precision

A&D Company, Ltd. (JAPAN)

URL: aandd.jp

A&D Engineering, Inc. (USA)

URL: andonline.com

A&D Australasia Pty Ltd. (Australia)

URL: andaustrolasia.com.au

A&D Instruments Ltd. (United Kingdom)

URL: andprecision.com

<German Sales Office>

URL: andprecision.com

A&D Korea Ltd. (South Korea)

URL: andk.co.kr

A&D Rus Co., Ltd. (Russia)

URL: and-rus.ru

A&D Instruments India (P) Ltd. (India)

URL: aanddindia.in

A&D Sciencetech Taiwan Ltd. (Taiwan)

URL: aandd.com.tw

A&D Instruments Thailand Ltd. (Thailand)

URL: thai.andprecision.com