

PA Tip Compatibility Table

Confirmed on March 5, 2021

Pipette	Tip maker	Model number	Compatibility	Remarks
PA-10/20	Gilson	D10	Compatible	
		DL10	Incompatible	Not possible to affix tip
	Eppendorf	epT.I.P.S. Standard for 10 uL (0030 000.811)	Compatible	
	Rainin	RC-10 / RCC-10	Compatible	
	Biohit	Optifit Tips for 0.1-10 uL	Compatible	
	Thermo	Finntip Flex 10	Compatible	
	BRAND	Pipette Tips for 0.1-20 uL	Compatible	
	Nichiryo	00-BMT-SS	Compatible	
	QSP	101	Compatible	
		114	Incompatible	Not possible to affix tip
	Watson	207C	Compatible	
	Nippon Genetics	37660	Compatible	
		10310	Compatible	
	Labcon	1038	Compatible	
AS ONE	V-10RE	Compatible		
PA-200	Gilson	D200	Compatible	
	Eppendorf	epT.I.P.S. Standard for 200 uL (0030 000.889)	Compatible	
	Rainin	RC-250 / RCC-250	Compatible	
	Biohit	Optifit Tips for 0.5-200 uL	Compatible	
	Thermo	Finntip Flex 200	Compatible	
	BRAND	Pipette Tips for 2-200 uL	Compatible	
	Nichiryo	00-BMT-SE	Compatible	It is necessary to adjust the height of the tip ejector
	QSP	110-NEW	Compatible	
		110-N	Compatible	
		110-B	Compatible	
	Watson	705Y	Compatible	
	Nippon Genetics	FG-301RS	Compatible	
		FG-302	Compatible	
		10340	Compatible	
Labcon	1093	Compatible		
AS ONE	V-200RNE	Compatible		
PA-1200	Gilson	D1000	Caution	At aspiration of 1200 uL, the tip holder becomes wetted Use is possible at aspiration of 1000uL
	Eppendorf	epT.I.P.S. Standard for 1000 uL (0030 000.927)	Compatible	
	Rainin	RC-1000 / RCC-1000	Compatible	
	Biohit	Optifit Tips for 10-1000 uL	Compatible	
	Thermo	Finntip Flex 1000	Compatible	
	BRAND	Pipette Tips for 50-1000 uL	Compatible	
	Nichiryo	00-BMT-L	Compatible	
		00-BMT2-LER	Compatible	
	QSP	111	Compatible	
	Watson	706B	Compatible	
	Nippon Genetics	34760	Compatible	
	Labcon	1045	Compatible	
	Greiner Bio-One	750251	Compatible	
	AS ONE	V-1000RNE	Compatible	
PA-10000	Gilson	D10 mL	Compatible	*1
	Eppendorf	epT.I.P.S. Standard for 10 mL (0030 000.765)	Compatible	
	Rainin	RC-10ML	Compatible	*1
	Thermo	Finntip 10 mL	Incompatible	Not possible to affix tip
	BRAND	Pipette Tips 1-10 mL	Compatible	*1
	Nichiryo	00-BMT2-Z	Incompatible	Not possible to affix tip
	QSP	097	Compatible	*1
	Watson	409C	Compatible	*1
Axygen	T-10ML-C	Compatible	*1	

\*1 MPA-10000 units produced before March 2019 cannot remove chips using the ejector.

**Compatibility**

Compatible	No problems with affixing tip, ejecting, or aspirating and dispensing The tip's accuracy has been confirmed to within +/-1.5% (at the maximum volume) of the levels of A&D's standard tips (With the MPA's volume calibration functionality, it is possible to increase accuracy further)
Caution	Limits to volumes for use
Incompatible	Not possible to affix tip, etc.

Compatibility has been confirmed by A&D, but it is not guaranteed in the future. Please refer to the compatibility table and check with your own chips.